

# Design high specificity CRISPR-Cas9 gRNAs: principles and tools

Heidi Huang, PhD






- 1** Introduction of CRISPR-Cas9
- 2** gRNA Design
- 3** Resources and Services
- 4** Q&A

# What is CRISPR?





- **CRISPR** – **C**lustered **R**egularly **I**nterspaced **S**hort **P**alindromic **R**epeats
- **Cas9** – **C**RISPR **a**ssociated **s**ystem. RNA-guided dsDNA-binding protein that has nuclease activity
- Type II CRISPR/Cas systems  
CRISPR locus: TracrRNA + Cas complex proteins + spacers/repeats



Small fragments of viral DNA are stored in CRISPR locus as **spacers**

Adapted from: Mali P. et al. **Cas9 as a versatile tool for engineering biology**. *Nat. Methods* (2013), 10(10):957-963


# What is CRISPR?



Adapted from: Mali P. et al. **Cas9 as a versatile tool for engineering biology.** *Nat. Methods* (2013), 10(10):957-963



# DNA Repair Enables Targeted Genome Editing





# Potential Applications for CRISPR-Cas9





## Genome editing

Wild-type Cas9





Cas9 Nickase




## Genome regulation, reorganization and visualization

Cas9<sub>nuclease-null</sub> Protein Fusions



Cas9<sub>nuclease-null</sub> Nucleic Acid

### Structural aggregation



Adapted from: Mali P. et al. **Cas9 as a versatile tool for engineering biology.** *Nat. Methods* (2013), 10(10):957-963

# Off-target of CRISPR-Cas9



Target site: 5'- NNNNNNNNNNNNNNNNNNNNNNNNNN **NGG**  
20 nt PAM


- ◆ CRISPR/Cas9 can tolerate 1-3 mismatches between gRNA and target sequence, which can lead to off target nuclease activity.
- ◆ Cas9 induces double stranded breaks, any off target nuclease activity can cause mutations in those genes, leading to possible oncogenesis.



# gRNA Design - The Goal



- ◆ Accurate gRNA target sites
- ◆ Lowest off-target potential
- ◆ Optimum target location



# gRNA Design – Step by Step



Step 1: Target gene analysis

Step 2: Find gRNA canonical sequences

Step 3: Off-target analysis

Step 4: Location analysis

Step 5: gRNA delivery

# Step 1: Target Gene Analysis



## ◆ Gene sequence analysis

- NCBI database
- Genome sequencing
- **Our recommendation: Always sequence the target gene in the cells/strains you work with before designing gRNA**

## ◆ Gene structure analysis

- Exons/CDS
- Introns

# Step 2: Find gRNA Canonical Sequences



## ◆ gRNA canonical sequences

5'- NNNNNNNNNNNNNNNNNNNNNN **NGG**  
20 nt PAM

## ◆ GN(20)GG

5'- GNNNNNNNNNNNNNNNNNNNN **NGG**  
20 nt PAM

## ◆ PAM sequences

- **NGG**
- NAG

# Step 3: Off-target Analysis



- ◆ Genome wide analysis to find similar sequences to GN(20)GG
  - Higher similarity, higher off-target risk

**Target site:** 5'- GTGTAACGGATAATGGAC **NGG**

Distal                      Seed                      PAM

- ◆ 3 criteria to gauge potential off-target sites
  - Existing of PAM sequence
    - NGG
    - NAG
  - Sequence similarity
    - higher similarity, higher risk
  - Sequence similarity at 3' increase off-target risk

# Step 4: Location Analysis



- ◆ Knock-out: NHEJ to induce indels
  - early exons
  - function domain
  
- ◆ Knock-in: HR with donor template
  - Both exon and intron can be targeted
  - CDS (Coding DNA sequence) on exon can be better than intron for donor design
  - **Avoid mis-targeting of donor template by gRNAs**

**Recommendation: 2 or more gRNAs for each target gene**

# Step 5: gRNA Delivery



## ◆ gRNA expression plasmid

- U6 promoter
- gRNA targeting region
- gRNA scaffold
- Termination signal



## ◆ RNA

## ◆ Transgenic

# Enhancing Specificity By Truncating gRNA Length



- ◆ gRNA sequences can be 17-20 nt in length to achieve similar levels of on-target gene editing
- ◆ Up to 10,000 fold improvement in target specificity when truncated (17 or 18 base pair) gRNA is used
- ◆ Using a shorter gRNA (17 or 18 nt) can greatly improve off-target specificity

Fu Y. et al. Improving CRISPR-Cas nuclease specificity using truncated guide RNAs. *Nat. Biotech.* (2014). 32:279-284




# Enhancing Specificity with Modified Nucleases



## ◆ Paired Cas9 nickases (Cas9n)


- D10A Cas9 mutant allows for single strand nicking



Adapted from Ran AF. et al. **Double nicking by RNA-guided CRISPR Cas9 for enhanced genome editing specificity.** *Cell* (2013). 154:1380-1389

## ◆ Dimeric CRISPR RNA-guided FokI nuclease (FokI-dCas9)

- FokI was fused to a catalytically inactive Cas9 (dCas9) mutant



Adapted from: Tsai SQ et al. **Dimeric CRISPR RNA-guided FokI nucleases for highly specific genome editing.** *Nat. Biotech.* (2014). 32:569-575

# Paired gRNA Design for Cas9n and FokI-dCas9



- ◆ Similar principle as gRNA design for wild-type Cas9
- ◆ Two gRNAs are needed to target each strand of the target gene respectively
- ◆ Off-set, the distance between two gRNAs, should be considered
  - Cas9n: -4 to 20bp, wider range -8 to 100bp
  - FokI-dCas9: 14-17bp

Ran AF. et al. **Double nicking by RNA-guided CRISPR Cas9 for enhanced genome editing specificity.** *Cell* (2013). 154:1380-1389  
Tsai SQ et al. **Dimeric CRISPR RNA-guided FokI nucleases for highly specific genome editing.** *Nat. Biotech.* (2014). 32:569-575

# gRNA Online Design Tool



## gRNA design tool for CRISPR genome editing

GenScript's bioinformatics scientists have developed a leading algorithm to design guide RNA (gRNA) for targeting Cas nucleases to your desired target gene. This free online gRNA design tool allows you to instantly determine the best gRNA sequences for your gene editing experiments. It's quick and easy to order expression-ready gRNA constructs. To get started, tell us about your target gene and the CRISPR/Cas system you plan to use:

**Option 1: Search by gene name.**

Enzyme:  ?

Human gene:  ?

Quickly compare and select from gRNA sequences designed to uniquely target your gene of interest. Our user-friendly display shows you where each gRNA aligns within the endogenous genome so you can choose which exons or specific loci to target. This tool leverages our genome-wide database of over 200,000 gRNA sequences pre-designed to meet stringent criteria for optimizing efficient and highly specific gene targeting.

**Option 2 (advanced): Search by sequence:**

Species:  ?

PAM:  ?

Enzyme:  ?

Gene sequence:

This tool performs a real-time calculation to identify Cas9-targetable loci within a FASTA sequence you enter, and shows potential off-target sites within the human reference genome. This tool is ideal for use to design gRNA targeting transgenes, mutations, or modified loci specific to certain cell lines whose genomes differ from the human reference genome. Caution: Because this tool does not require that your input sequence be found in the human reference genome, it will not recognize exon boundaries or errors in your input.

- ◆ WT-Cas9 and Cas9n
- ◆ Multiple input options
- ◆ Human, mouse, CHO

<http://www.genscript.com/gRNA-design-tool.html>

# gRNA Online Design Tool



## TP53 CRISPR guide RNA, tumor protein p53 CRISPR guide RNA[human]

gRNA (name)	Exon	gRNA sequence	On-target sequence	Highest Risk Score ?	Potential off-target site with highest risk score ?	Price (per gRNA sequence)	Select
TP53 CRISPR guide RNA 1	1	CCACCGTC CAGGGAGC AGGTAGC	NC_000017.11 7687475..7687453	0	not found	\$199.00	<input type="checkbox"/>
			GCTACCTGCTCCCTGGACGGTGG       CCACCGTCCAGGGAGCAGGTAGC				
▪ ▪ ▪							
TP53 CRISPR guide RNA 4	1	GGACACTTT GCGTTTCGG GCTGGG	NC_000017.11 7687439..7687417	0	not found	\$199.00	<input type="checkbox"/>
TP53 CRISPR guide RNA 5	10	GCATGGGC GGCATGAAC CGGAGG	NC_000017.11 7674238..7674216	0	NC_000003.11 168368424..168368446 GCATGGGCGGCATGAACCGGAGG       TGATGGGCAGCATGAACCGGATG	\$199.00	<input type="checkbox"/>
			GCATGGGCGGCATGAACCGGAGG       GCATGGGCGGCATGAACCGGAGG				
TP53 CRISPR guide RNA 6	6	GGGCAGCT ACGGTTTCC GTCTGG	NC_000017.11 7676057..7676035	47	NC_000005.9 179022511..179022533 GGGCAGCTACGGTTTCCGTTCTGG       GGGCAGCTACGGTTTCTTTTTGG	\$199.00	<input type="checkbox"/>
			GGGCAGCTACGGTTTCCGTTCTGG       GGGCAGCTACGGTTTCCGTTCTGG				
TP53 CRISPR guide RNA 7	6	GGCAGCTA CGGTTTCC GTCTGGG	NC_000017.11 7676056..7676034	53	NC_000002.11 203053552..203053574 GGCAGCTACGGTTTCCGTTCTGGG       CCAAAGCTACGGTTTCTGTCTTGG	\$199.00	<input type="checkbox"/>
			GGCAGCTACGGTTTCCGTTCTGGG       GGCAGCTACGGTTTCCGTTCTGGG				


- ◆ Off-target risk scores (0- 58) by analyzing genome-wide sequence similarity
- ◆ Risk score <49 is recommended.

# GenCRISPR™ gRNA Construct Service



## GenCRISPR™ gRNA constructs

- free gRNA design
- ready-to-use plasmids for **\$199**



	Service Steps	Deliverables	Timeline	Price
GenCRISPR gRNA construct service	<ul style="list-style-type: none"><li>• gRNA design</li><li>• gRNA synthesis</li><li>• Cloning into vector</li><li>• Validation by PCR and enzyme digestion</li><li>• Validation by sequencing</li></ul>	<ul style="list-style-type: none"><li>• 4ug plasmid for each gRNA construct</li><li>• Final report with QC data</li></ul>	10 business days	<b>\$199/construct</b>

- ◆ Full service covers both design and construction
- ◆ \$199/construct, lowest price on the market
- ◆ Free standard gRNA vector and complementary human codon-optimized Cas9 expression vector
- ◆ **All-in-one vectors are available now!**

# GenCRISPR™ Cell Line Service



## GenCRISPR™ Genome Editing Service

- **30% off** from now until November 30<sup>th</sup> 2014!



## Service details

GenCRISPR™ Cell line Service	Deliverables	Timeline	Pricing
<b>Knock-out cell line, single clone (SC1652)</b> Customer specifies gene/locus region and cell line of interest*	<ul style="list-style-type: none"> <li>• Single clone, target sequence validated</li> <li>• Detailed report with gRNA targeting sequence design and sequencing data of target region</li> </ul>	<b>9-10 weeks</b> , depending on complexity of target and cell line growth	<b>\$8000</b>
<b>Knock-in cell line, single clone (SC1663)</b> Customer specifies target gene/locus and provides insert sequence or mutation on target Customer also specifies cell line of interest*	<ul style="list-style-type: none"> <li>• Single clone, target sequence and knock-in gene sequence validated</li> <li>• Detailed report with gRNA targeting sequence design, and knock-in homology arm sequence</li> </ul>	<b>12 weeks</b> , depending on complexity of target and insertion sequence, and cell line growth	<b>\$9000</b>

# Summary



- ◆ CRISPR-Cas9 is an efficient and easy to implement system for genome editing.
- ◆ Cas9 can tolerate 1-3 mismatches and generate off-target mutations.
- ◆ gRNA design steps: Target gene analysis, find canonical gRNA sequences, off-target analysis, location analysis, and gRNA delivery.



# Summary



- ◆ Modified Cas9 nucleases, including Cas9n and FokI-dCas9, improve targeting specificity. Paired gRNAs are needed.
- ◆ GenScript offers free gRNA online design tool, featuring off-target risk ranking.
- ◆ GenScript offers GenCRISPR™: a complete genome editing solution including gRNA design and construction, and custom cell line development.





*Thank you for your participation*  
*We wish you all success in your research*  
**Email me: [Heidi.Huang@GenScript.com](mailto:Heidi.Huang@GenScript.com)**

Register for other webinars in the GenScript Webinar Series @ <http://www.genscript.com/webinars.html>



**November 12, 2014/8:00 am, 2:00pm EST**

**Identify the optimal protein purification strategy for your recombinant protein production**



**November 20, 2014/8:00 am, 2:00pm EST**

**Anti-idiotypic antibodies – A powerful tool for antibody drug studies**