

Pfu DNA Polymerase

Cat. No.: E00021

Size: 250 U

Unit Definition:

1 unit of enzyme incorporates 10 nmol of dNTP into acid-insoluble material in 30 min at 74 °C.

Version: 08/03/2006

PRODUCT INFORMATION

Description

Pfu DNA Polymerase is a thermostable DNA Polymerase isolated from an *E.coli* strain that carries the pol gene from *Pyrococcus furiosus*. As one of the most thermostable DNA Polymerase known with 3' to 5' proofreading activity, the enzyme catalyze DNA synthesis at optimal temperature near 75°C with very low error rate (about eight times more accurate than *Taq* DNA Polymerase). *Pfu* DNA Polymerase generates blunt-ended PCR fragments that are required for blunt-end PCR cloning.

10 X reaction Buffer (with Mg²⁺)

200 mM Tris-HCl (pH 8.7 at 25°C), 100 mM (NH₄)₂SO₄, 100 mM KCl, 1% (v/v) Triton X-100, 1 mg/ml BSA and 20 mM MgSO₄.

Storage Buffer and Concentration

Supplied in 2.5 units/μl in 50 mM Tris HCl (pH 8.1), 0.1 mM EDTA, 1 mM DTT, 0.1% (v/v) Nonidet P40, 0.1% (v/v) Tween 20 and 50% (v/v) glycerol.

Storage

-20 °C.

Applications

Pfu DNA Polymerase can be used in most applications including the following:

- High fidelity PCR amplification.
- High fidelity blunt-end cloning.
- Site-directed mutagenesis.

QC Tests

PCR⁺ performance and activity.




Figure 1. PCR amplification of a 544 bp gene from *E. coli* genomic DNA using *Pfu* DNA Polymerase. 0.25, 0.50, 0.75, 1.00 and 1.25 U of the enzyme were used in 50 μl of PCR reaction 1 (Lane 1), 2 (Lane 2), 3 (Lane 3), 4 (Lane 4) and 5 (Lane 5), respectively. 5 μl of PCR reaction was loaded in each lane. M is a 100 bp DNA Marker.




Figure 2. GenScript *Pfu* was compared with the *Pfu* from a leading brand (A) for the PCR amplification of a 544 bp gene from *E. coli* genomic DNA. 0.4, 0.6, 0.8, and 1.0 U of the enzymes were used in 50 μl of PCR reactions as shown in Figure 2. 10 μl of PCR reaction was loaded in each lane. M is a 100 bp DNA Marker.

** For non-clinical research use only. **

Tel: 732-885-9188 Fax: 732-210-0262 E-Mail: info@genscript.com www.genscript.com