

# LApró Taq DNA Polymerase



Technical Manual No. TM0219

Version 01152007

I	Description.....	1
II	Applications.....	1
III	Contents .....	1
IV	Activity .....	1
V	Unit Definition.....	2
VI	Storage Buffer and Concentration.....	2
VII	Storage.....	2
VIII	General Protocol Using LApro Taq DNA Polymerase .....	2
IX	Troubleshooting .....	3
X	Order Information.....	4

## I. DESCRIPTION

LApró Taq DNA Polymerase is a polymerase possessing 5'-3' DNA polymerase activity that is optimized to amplify longer GC-rich fragments and 3'-5' proof-reading activity. As one of the most thermostable DNA polymerase known with 3'-5' proofreading activity, the enzyme catalyze DNA synthesis at optimal temperature near 75°C with very low error rate (about ten times more accurate than standard Taq DNA Polymerase). The enzyme is optimized to achieve amplification of 17.5 kb products from human genomic DNA and 20 kb from λDNA. LApro Taq DNA Polymerase in long PCR\* amplifications can reduce smearing and virtually eliminate unwanted background artifacts. LApro Taq DNA Polymerase generates A-ended PCR fragments that are required for TA cloning.

## II. APPLICATIONS

LApró Taq DNA Polymerase can be used in most applications including the following:

- Up to 17.5 kb from human genomic DNA
- High fidelity amplification
- More stable than normal Taq DNA Polymerase
- TA cloning
- Site-directed mutagenesis
- Fragment analysis

## III. CONTENTS

Components	E00013	E00014
LApró Taq DNA Polymerase	250 U	250 U
2x PCR Buffer	2X1.5 ml	2X1.5 ml
dNTP, 10 mM (D0056)		50ul
Control template & primers		25 rxn

## IV. ACTIVITY

LApró Taq DNA Polymerase shows very good fidelity and an average error rate around  $2 \times 10^{-6}$ . LApro Taq DNA Polymerase has 5'-3' exonuclease activity and 3'-5' proofreading activity. It also has the extendase activity allowing TA cloning.

## V. UNIT DEFINITION



## Lapro Taq DNA Polymerase

One unit of enzyme incorporates 10 nmol of dNTP into acid-insoluble material in 30 minutes at 74°C.

### VI. STORAGE BUFFER AND CONCENTRATION

Supplied in 5 units/μl in 50 mM Tris-HCl (pH8.1), 0.1 mM EDTA, 1 mM DTT, 0.1% (v/v) Nonidet P40, 0.1% (v/v) Tween 20, and 50% (v/v) glycerol.

### VII. STORAGE

Store at -20°C. Stable for at least one year, if stored properly.

### VIII. GENERAL PCR PROTOCOL USING LAPRO TAQ DNA POLYMERASE

This is a general PCR amplification protocol, optimization may be needed to get satisfactory results.

1. Set up 50 μl PCR reaction by adding the following reagents to a thin-walled PCR microcentrifuge tube or plate and mixing gently.

Reagent	Volume	Final Concentration
Water, PCR grade	20 μl	
2x Buffer	25 μl	1x
10 mM dNTP	1 μl	200 nM each
20 μM Forward primer	1 μl	400 nM
20 μM Reverse primer	1 μl	400 nM
Template DNA	1 μl	1-100 ng/reaction
5 U/μl LApr Taq DNA Polymerase	1 μl	0.1 U/μl
Total	50 μl	

Note: when use our positive control in complete system (E00014) just add 50ul ddH<sub>2</sub>O to the control tube and fully dissolved (gently), then tip 2ul to the PCR microcentrifuge tube as template and primers. Pcr product length is 20kb.

2. Suggested PCR profile:
  - 3 Step PCR
    - Initial Denaturation: 94-96°C for 2 min.
    - 25-30 Cycles of Denaturation: 94-96°C for 30 sec.
    - Annealing: 50-60°C for 30 sec. (primer T<sub>m</sub> - 5°C)
    - Extension: 72°C for 30 sec to 5 min. (1.5-2.5 kb/min)
    - Final Extension: 72°C for 7 min.
  - 2 Step PCR
    - Initial Denaturation: 94-96°C for 1 min.
    - 20-30 Cycles of Denaturation: 94-96°C for 20 sec.
    - Annealing & Extension: 68°C for Xmin. (1-2.5 kb/min)
    - Final Extension: 72°C for 7-15 min.
3. Examples of Using LApr Taq DNA Polymerase



## Lapro *Taq* DNA Polymerase




Figure 1. PCR amplification of 1 kb~20 kb fragments from  $\lambda$ DNA using LAprro *Taq* DNA Polymerase. 1U of enzyme was used in 50  $\mu$ l of PCR reaction. Lane 1 - 1kb, Lane 2 - 2kb, Lane 3 - 3kb, Lane 4 - 4kb, Lane 5 - 5kb, Lane 6 - 6kb, Lane 7 - 8kb, Lane 8 - 10kb, Lane 9 - 12kb, Lane 10 - 20kb. 5 $\mu$ l of PCR reaction was loaded in each lane. M is DNA KB Ladder (GenScript, Cat. No. M101R).




Figure 2. PCR amplification of 17.6kb fragments (lane 2, 3, 4) from human genomic DNA and 20kb (lane 6, 7, 8) from  $\lambda$ DNA using LAprro *Taq* DNA Polymerase.

Lane 1: 5 U standard *Taq* DNA Polymerase were used to amplify 17.6kb fragment from human genomic DNA

Lane 2, 3 & 4: 1 U, 2 U and 5 U of LAprro *Taq* DNA Polymerase were used in 50  $\mu$ l of PCR reaction to amplify 17.6 kb fragment from human genomic DNA

Lane 5: 5 U KOD DNA Polymerase were used to amplify 20 kb fragment from  $\lambda$ DNA

Lane 6, 7 & 8: 1 U, 2 U and 5 U of LAprro *Taq* DNA Polymerase were used in 50  $\mu$ l of PCR reaction to amplify 20 kb fragment from  $\lambda$ DNA. 3  $\mu$ l of PCR reaction was loaded in each lane. M is a  $\lambda$  *Hind*III DNA Marker.

## IX. TROUBLESHOOTING

The table below is guideline for troubleshooting.

Problem	Probable Cause	Solution
No PCR DNA	Extension time is too short; One or more PCR components	LAprro <i>Taq</i> DNA Polymerase requires a minimum extension time of 1 minute for 1~2 kb template. Run a positive control to make sure the PCR



## Lapro Taq DNA Polymerase

	<p>may be missing; PCR conditions are not optimized; the annealing temperature may be too high; more cycles may be needed; the denaturation time may be too short;</p> <p>The primers may not be designed optimally.</p> <p>Target template is highly GC-rich.</p>	<p>reaction is properly set-up. Optimize the PCR conditions by decreasing annealing temperature in 2-4°C increments, or increasing the number of cycles, or increasing the denaturation time in 10 second increments. It is recommended to change only one parameter each time.</p> <p>The primer designing is critical for high quality PCR. Longer primers of 25-30 nucleotides with a GC content of 45-60% and with a more stable 5'-end than 3'-end usually make good primers.</p> <p>The target will be difficult to denature even with a longer denaturation step.</p>
Non-specific DNA products	<p>The primers may not be designed optimally. Annealing temperature is too low.</p>	<p>See above. Optimize the PCR conditions by increasing annealing temperature in 2-4°C increments, or decreasing the number of cycles.</p>
False positive	<p>Reagents are contaminated.</p>	<p>It is recommended that a negative control without using genomic DNA be run to make sure no contamination occurs.</p>

### X. ORDER INFORMATION

Lapro Taq DNA Polymerase, Cat. No. E00013

Lapro Taq DNA Polymerase Complete System, Cat. No. E00014

Telephone: 732-885-9188, 732-357-3839

Fax: 732-210-0262, 732-885-5878

E-mail: [info@genscript.com](mailto:info@genscript.com)

### For Research Use Only.

#### GenScript Corporation

120 Centennial Ave., Piscataway, NJ 08854

Tel: 732-885-9188, 732-357-3839

Fax: 732-210-0262, 732-885-5878

E-mail: [info@genscript.com](mailto:info@genscript.com)

Web: <http://www.Genscript.com>

\* The PCR process is covered by US. Patent numbers 4683195 and 4683202 issued to Cetus and owned by Hoffman-La Roche Inc. GenScript does not encourage or support the unauthorized use of the PCR process. Use of this product is recommended for persons that either have a license to perform PCR or are not required to obtain a license.