

I	Description.....	1
II	Key Features.....	1
III	Contents.....	1
IV	Stability.....	2
V	Storage.....	2
VI	Storage.....	2
VII	General PCR Protocol Using 2X Taq DNA Master Mix.....	2
VII	Order Information.....	3

## I. DESCRIPTION

2X Taq Master Mix is a premixed 2X concentrated solution of Taq DNA Polymerase (GenScript, Cat. No. E00007), reaction buffer, MgCl<sub>2</sub> and dNTPs. 2X Taq Master Mix contains all components for PCR\*, except DNA template and primers. The mixture is optimized for consistent and efficient routine PCR amplifications. It can amplify up to 8 kb fragment from lambda DNA. For a 50 µl reaction, simply add 25 µl of 2X Taq Master Mix to primers, DNA template and PCR-Qualified H<sub>2</sub>O.

## II. KEY FEATURES

- Taq DNA Polymerase in ready-to-use mix
- Low contamination risk
- Low risk of pipetting errors

## III. CONTENTS

- 0.1 U/µl Taq DNA Polymerase (GenScript, Cat. No. E00007)
- Reaction buffer
- 3 mM MgCl<sub>2</sub>
- 0.4 mM dNTPs

## IV. STABILITY

- Freeze-thaw stability of 2X *Taq* master mix: Following 25 freeze-thaw cycles, no effect on performance is observed.




Fig. 1 Stability after 25 freeze-thaw cycles.

2 kb fragment amplification.

- Lane 1 25  $\mu$ l of 2X *Taq* Master Mix
- Lane 2 25  $\mu$ l of 2X *Taq* Master Mix after 25 freeze-thaw cycles
- Lane 3 2.5 U *Taq* DNA Polymerase (GenScript, Cat. No. E00007)

- Stability at 4°C: No effect on performance is observed after storage at 4°C for 2 months.




Fig. 2 Stability at 4°C

0.5 kb fragment amplification.

- Lane 1 25  $\mu$ l of 2X *Taq* Master Mix
- Lane 2 25  $\mu$ l of 2X *Taq* Master Mix storage at 4°C for 2 months
- Lane 3 2.5 U *Taq* DNA Polymerase (GenScript, Cat. No. E00007)

## V. STORAGE

This product is shipped on blue ice. Store the product at -20°C.

## VI. GENERAL PCR PROTOCOL USING 2X TAQ MASTER MIX

This is a general PCR amplification protocol, optimization may be needed to get satisfactory results.

- Thaw the 2X *Taq* Master Mix at room temperature. Vortex the 2X *Taq* Master Mix and then spin it briefly in a microcentrifuge to collect the material in the bottom of the tube.
- Prepare one of the following reaction mixes on ice:

- **For a 25  $\mu$ l reaction volume:**

Component	Volume	Final Concentration
2X <i>Taq</i> Master Mix	12.5 $\mu$ l	1X
Upstream Primer, 10 $\mu$ M	0.5 $\mu$ l	0.1–1.0 $\mu$ M
Downstream Primer, 10 $\mu$ M	0.5 $\mu$ l	0.1–1.0 $\mu$ M
DNA Template	1-5 $\mu$ l	<500 ng
Nuclease-Free Water to	25 $\mu$ l	

- **For a 50  $\mu$ l reaction volume:**

Component	Volume	Final Concentration
2X <i>Taq</i> Master Mix	25 µl	1X
Upstream Primer, 10 µM	1 µl	0.1–1.0 µM
Downstream Primer, 10 µM	1 µl	0.1–1.0 µM
DNA Template	1-5 µl	<500 ng
Nuclease-Free Water to	50 µl	

- Gently mix the reaction and spin down in microcentrifuge.
- Set up cycling conditions for a routine PCR reactions:
  - Initial Denaturation 94-95°C for 1-5 minutes
  - 25-40 cycles 94-96°C for 30 sec.  
45-70°C for 10-30 seconds  
72°C for X minutes (1 min/kb)
  - Final extension: 72°C for 7 minutes
  - Final soak 4-10°C

## VII. ORDER INFORMATION

2X *Taq* Master Mix, Cat. No. E00019

**For Research Use Only.**

### GenScript USA Inc.

860 Centennial Ave., Piscataway, NJ 08854

Tel: 732-885-9188, 732-885-9688

Fax: 732-210-0262, 732-885-5878

E-mail: [info@genscript.com](mailto:info@genscript.com)

Web: <http://www.Genscript.com>

\* The PCR process is covered by US. Patent numbers 4683195 and 4683202 issued to Cetus and owned by Hoffman-La Roche Inc. GenScript does not encourage or support the unauthorized use of the PCR process. Use of this product is recommended for persons that either have a license to perform PCR or are not required to obtain a license.