



I	Description	1
II	Key Features	1
III	Kit Contents	1
IV	Storage.....	2
V	Protocol	2
VI	Related Products.....	3
VII	Ordering Information.....	3

I. DESCRIPTION

The **Caspase-8 Colorimetric Assay Kit (Cat. No. L00303)** provides a simple and convenient means of assaying the activity of caspase-8. The activation of ICE-family proteases and caspases initiates apoptosis in mammalian cells. The assay for caspase-8 is based on spectrophotometric detection of the chromophore *p*nitroanilide (*p*NA) after cleavage from the labeled substrate IETD-*p*NA. The *p*NA light emission can be quantified using a spectrophotometer or a microtiter plate reader at 400 or 405 nm. By comparing the absorbance of *p*NA from an apoptotic sample to an uninduced control, the relative increase of caspase-8 activity can be determined.

II. KEY FEATURES

- **Simplified Procedure:** The components in this kit are provided ready-to-use.
- **Supreme Reliability:** The kit produces highly reproducible results.

III. KIT CONTENTS

The Caspase-8 Colorimetric Assay Kit is available (L00303) employs a substrate IETD-*p*NA. Each kit contains enough reagents for one hundred apoptosis assays.

Components	Cat. No. L00303 100 assays
Lysis Buffer	15.0 ml
2 Reaction Buffer	5.0 ml
Caspase-8 Substrate	500 µl, store away from light
DTT	150 µl

IV. STORAGE

Store the Lysis Buffer and 2X Reaction Buffer at 4°C after opening. Store the rest of the kit at -20°C. All reagents will remain stable for one year.

V. CASPASE-8 COLORIMETRIC ASSAY KIT PROTOCOL

A. General Considerations

- Aliquot enough 2X Reaction Buffer for the number of assays to be performed. Add DTT into the 2X Reaction Buffer immediately to a final concentration of 10 mM before use (10 µl of 1.0 M DTT stock in per 1 ml 2X Reaction Buffer).
- After thawing, store the Lysis Buffer at 4°C.
- Do not expose the Caspase-8 Substrate to light.

B. Assay Procedure

1. Induce apoptosis and set a negative control at the same time.
2. Control at the same time.
3. Collect the cells and wash them twice with PBS by centrifugation at 2000 rpm for five minutes.
4. Collect 3×10^6 - 5×10^6 cells and try to remove the PBS.
5. Add 50 µl cold prepared Lysis Buffer into the collected cells and mix well. (Note: Before using the Lysis Buffer, add 0.5 µl DTT to 50 µl Lysis Buffer.)
6. Incubate the cells on ice for 20-60 minutes, and vortex with vibration three or four times for 10 seconds each time or freeze and thaw the cells two or three times.
7. Centrifuge at 10,000 rpm at 4°C for one minute.
8. Remove the supernatant and transfer it to a new tube. Put the tube on ice.
9. Obtain a small quantity of supernatant to assay the protein concentration by the Bradford or BCA method.
10. Take 50 µl supernatant containing 100-200 µg protein. If there is not enough volume, add Lysis Buffer.
11. Add 50 µl of 2X Reaction Buffer to each sample. (Note: Add 0.5 µl DTT to 50 µl 2X Reaction Buffer.)
12. Add 5 µl Caspase-8 Substrate (200 µM final concentration) and incubate at 37°C for four hours, away from light.
13. Determine the extinction value of the samples using a spectrophotometer (100 µl micro quartz cuvet) or a microtiter plate reader at 400 nm or 405 nm. Obtain the results of the induced group's caspase-8 activity by computing $OD_{\text{inducer}}/OD_{\text{negative control}}$.

Note: The background from cell lysates and buffers from the readings of induced samples and negative controls should be subtracted before calculating the increase in caspase-8 activity.

VI. RELATED PRODUCTS**Caspase-3 Colorimetric Assay Kit, Cat. No. L00289****Caspase-9 Colorimetric Assay Kit, Cat. No. L00304****Caspase-2 Colorimetric Assay Kit, Cat. No. L00294****Caspase-6 Colorimetric Assay Kit, Cat. No. L00302****VII. ORDERING INFORMATION****Caspase-8 Colorimetric Assay Kit, Cat. No. L00303**

GenScript USA Inc
860 Centennial Ave., Piscataway, NJ 08854
Toll-free Tel: 1-877-436-7274
Fax: 1-732-210-0262
E-mail: product@genscript.com
Web: www.genscript.com

For Research Use Only.