

## PRODUCT INFORMATION

Version: 07/29/2008

### Protein Marker for Fluorescent Western Blotting

**Cat. No.** M00124

**Size:** 100 Lanes

**Category:** Protein Maker

#### Description

This protein marker is designed for convenient protein band identification in fluorescent Western Blotting. It is a lyophilized mixture of five recombinant proteins that are able to bind to primary and/or secondary antibodies derived from a wide range of host species. It thus enables direct visualization of both the protein marker and users' samples on the same western blot membrane without any additional reagents.

#### Calculations

The apparent molecular weights of its five protein bands are 22 kDa, 40 kDa, 60 kDa, 85 kDa and 120 kDa.

#### Reconstitute

Reconstitute the powder with 500  $\mu$ l reconstitution buffer.

Load 5  $\mu$ l per lane for a mini-gel.


#### Storage

After reconstitution, aliquot and store at  $-20^{\circ}\text{C}$  or below for long-term storage. Aliquots may be stable at  $4^{\circ}\text{C}$  for one month.

#### Notes:

1. Non-BSA based blocking reagents (such as 5% skimmed milk) are recommended for blocking.
2. For best results, the loading volume of the marker should be optimized for different experiment conditions (from 2.5  $\mu$ l to 10  $\mu$ l).
3. This marker binds to the constant region of IgG, so protein A or protein G fluorophore conjugates cannot be used as the secondary antibody.

#### Example



10  $\mu$ l and 5  $\mu$ l of reconstituted protein marker for Fluorescent Western Blotting were loaded in 12% Express SDS-PAGE gel (GenScript, MG012W10). Antibodies used were Rabbit total IgG and IRDye 680 conjugated Goat Anti-Rabbit IgG (LI-COR) (Fig.A) and Mouse total IgG and IRDye 800CW conjugated Goat Anti-Mouse IgG (LI-COR) (Fig.B). Images were acquired using Odyssey Infrared Imaging system (LI-COR).