



I	Blocking	1
II	Incubation with Primary Antibody	1
III	Incubation with Second Antibody	1
IV	Chemiluminescent Reaction	1
V	Peptide Competition	1
VI	Ordering information	2

General Protocol for Western Immunoblotting

I. BLOCKING

1. Block membrane by incubating 1 hour at room temperature or overnight at 4°C with shaking in Blocking Solution (5% BSA, 0.05% Tween-20 in TBS (50mM Tris, 100mM NaCl, pH 7.6).

II. Incubation with Primary Antibody

2. Dilute primary antibody at the appropriate dilution in Blocking Solution.
 3. Incubate the membrane with diluted primary antibody for 1 hour at 37°C, or 2 hours at room temperature, or overnight at 4°C with agitation.
 4. Remove antibody solution. Wash the membrane 3 times for 5-10 minutes each time at room temperature in TBST (50mM Tris, 100mM NaCl, 0.05% Tween-20, pH 7.6) with shaking. Note: Increase the concentration of Tween-20 to 0.1% reduces the background and increases the specificity, but it will reduce the sensitivity.

III Incubation with Second Antibody

5. Incubate membrane with secondary HRP conjugate diluted (according to manufacturer's instructions) in Blocking Solution for 1 hour at room temperature with shaking.
 6. Repeat Step 4.
 7. Wash membrane with TBS for 2-5 minutes before proceeding Chemiluminescent Reaction.

IV. Chemiluminescent Reaction

8. Prepare and use the Chemiluminescent substrate according to the manufacturer's instructions.
 9. Immediately wrap the membrane and expose to X-ray films for 10 second to 1 hour period. The exposure time may vary according to the amount of antibody and antigen.

V Peptide Competition

Before proceeding Western Immunoblotting, add Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.



VI. ORDERING INFORMATION

Telephone: 732-885-9188

Fax: 732-210-0262, 732-885-5878

Email: info@genscript.com

For Research Use Only.

GenScript Corporation

120 Centennial Ave., Piscataway, NJ 08854

Tel: 732-885-9188

Fax: 732-210-0262, 732-885-5878

Email: info@genscript.com

Web: www.genscript.com