

# Anti-Idiotypic Antibodies: Powerful Tools for Antibody Drug Development

Michelle Parker, Ph.D.

[Michelle.parker@genscript.com](mailto:Michelle.parker@genscript.com)



# Table of Contents




- 1 What is an Anti-Idiotype Antibody?
- 2 Anti-Idiotype Antibody Applications
- 3 Obstacles & Solutions to the Generation of Anti-Idiotype Abs
- 4 Downstream Assay Development
- 5 Features of GenScript's Anti-Idiotype Antibody Services
- 6 GenScript Anti-Idiotype Antibody Packages

# Antibody: Structure and Function



**Antibody (Ab):** *Recognition proteins found in the serum and other bodily fluids of vertebrates that react specifically with the antigens that induced their formation.*


- ◆ Overall structure:
  - 2 identical light chains (blue)
  - 2 identical heavy chains (green/purple)
- ◆ Variable regions and constant regions
- ◆ 5 classes of Abs:
  - IgG, IgA, IgM, IgD, IgE
  - All contain either  $\lambda$  or  $\kappa$  light chains
  - Biological effector functions are mediated by the C domain
- ◆ Chemical structure explains 3 functions of Abs:
  1. Binding versatility
  2. Binding specificity
  3. Biological activity




# Antibody Binding Regions




- ◆ **Idiotope** – the antigenic determinants in or close to the variable portion of an antibody (Ab)
- ◆ **Paratope** – the part of an Ab that recognizes an antigen, the antigen-binding site of an Ab or **complementarity determining region (CDR)**
- ◆ **Epitope** – the part of the antigen to which the paratope binds




# Anti-Idiotypic Antibodies



- ◆ Anti-idiotypic antibodies (Anti-IDs) – Abs directed against the paratope (or CDR region) of another Ab



Anti-Idiotypic  
Directed against CDR regions




Hypervariable regions (or the idiotype of an Ab)  
are the antigenic determinants

# Different types of anti-IDs




## Complex-specific




- Drug-target complex specific
- Not inhibitory
- Detects bound drug only

## Antigen-blocking




- Paratope-specific
- Inhibitory
- Neutralizing
- Detects free drug

## Non-blocking



- Not paratope-specific
- Not inhibitory
- Detects total drug (free, partially bound, fully bound)



## Why use anti-IDs?

- ◆ Ideal for pre-clinical research and antibody drug development studies
- ◆ Allow monitoring of therapeutic Abs in samples
- ◆ Allow detection of Ab biotherapeutics that closely resemble circulating human Ig



**Anti-IDs, specific for the unique variable region of the therapeutic Ab, are ideal for this purpose**




- ◆ Anti-IDs can be powerful reagents for the following applications:
  1. **Pharmacokinetic (PK) studies**
    - Used to measure the drug level in patient samples
  2. **Immunogenicity (anti-drug antibody) assays**
    - Used as a positive control or reference standard





## 3. Vaccine Development




General flow chart of anti-idiotypic antibody network




# Antibody Production Overview

# How are Monoclonal Antibodies Made?



# How are Polyclonal Antibodies Made?



# Challenges During Anti-ID Production




## 1) CDR region of Ab may not be very immunogenic

- **Antigen immunogenicity:** *The ability of a particular molecule to elicit an immune response determined by whether the immune system can recognize the antigen*

## 2) Percentage in antiserum very low

- Most Abs raised target the Fc region
- Anti-ID clone may be missed if conventional screening methods are used

## 3) Must optimize clone selection from numerous anti-IDs






## Overcoming the Obstacles to Anti-ID Production

# 1. CDR Region Not Immunogenic



## Solution

- ◆ Use antibody drug-KLH conjugate as immunogen
- ◆ Use  $F(ab)_2$  (fragment, antigen binding) as immunogen
  - Constant heavy chain determinants are absent
- ◆ Adjust immunization schedule
- ◆ Consider different animal host



## 2. Low Percentage of Anti-IDs in Serum



### Solution

- ◆ Use proper secondary Ab (anti-mouse Fc to prevent cross-reaction with human Ig)
- ◆ HTP binding screening
- ◆ Capture ELISA to select natural epitope recognized by Ab






# 3. Optimize Clone Selection from Numerous Anti-IDs



## Solution

- ◆ Antibody pairing  
(for PK study)
- ◆ Affinity ranking  
( $K_{\text{off}}$  and  $K_{\text{d}}$  ranking)
- ◆ Epitope binning  
(ELISA or SPR-based)
- ◆ Antigen blocking


## Antigen blocking assay



# Optimizing clone selection: Epitope binning



- ◆ Abs tested in a pairwise combinatorial manner
- ◆ Abs that compete for binding to the same binding region are grouped into the same epitope bin.




Mouse anti-ID  
(against epitope 2)

Mouse anti-ID-HRP  
(against epitope 1)

Target human Ab  
coated plate

Same epitope results in ↓ signal



## Epitope binning by Octet RED 96



## Case Study:

Assess relative affinity rank	Epitope			
	Epitope 1	Epitope 2	Epitope 3	Epitope 4
1	4F5	5A7	12A7	9A2
2		5B11	3E5	19C12
3		3H9		
4		3B3		
5		5G7		

# GenScript Anti-ID Development Protocol






## Downstream Assay Development

# Anti-IDs Offer Flexibility for Specialized Assay Development




- ◆ Due to different binding modes and levels of affinity
- ◆ Using a combination of anti-IDs one can obtain:
  1. PK profile
  2. Degree of target saturation by the antibody drug
- ◆ Possible binding modes of anti-IDs:
  - **Complex-specific** anti-IDs detect **bound** Ab drugs directly
  - **Antigen-blocking** (neutralizing) anti-IDs detect **free** drug
  - **Non-blocking** (non-inhibitory) anti-IDs detect **free and bound** drug


Antigen-blocking



Non-blocking



Complex-specific



# Guide for how to use concentration data



Concentration data	Intended use
$mAb_{\text{free}}$ or $mAb_{\text{total}}$	To characterize in vivo PK behavior and to project human PK To assess PK and PD response relationships in animal models To calculate the safety margin and determine the safe starting dose or efficacious dose
$L_{\text{free}}$	To assess the inhibition of $L_{\text{free}}$ following drug treatment; understand the dynamics of $L_{\text{free}}$
$L_{\text{total}}$	To assess the redistribution/modulation of $L_{\text{total}}$ following drug treatment; understand the dynamics of $L_{\text{total}}$
$L_{\text{free}}/L_{\text{total}}$ ratio	To compare in vivo and in vitro binding affinities
$mAb$ and $L$ dynamic relationship	To understand the underlying dynamic relationship between $mAb$ and $L$ to facilitate dose and dosing schedule selection

**L** = target ligand in circulation

**mAb** = monoclonal antibody therapeutic

Table adapted from The AAPS Journal, Vol.13 (1), 99-110, March 2011.



1

**Pharmacokinetic Assays**

2

**Immune Response Assays**

# Pharmacokinetic (PK) Assays




- ◆ Purpose: Detect and quantitate human Ab drugs in serum
- ◆ Human serum contains **5-9.5 mg/ml IgG1**
- ◆ Required sensitivity for Ab drug is in the **ng/ml range**
- ◆ Thus, assays must accurately detect Ab drug despite a million-fold excess of similar molecules
- ◆ Possible agents for detection include:
  - Antigen (i.e. drug target)
  - **Anti-ID antibodies**




# Typical PK Assay Formats




## 1 Antigen Capture




## 2 Anti-ID-Bridging



## 3 Anti-ID Capture Sandwich



## 4 Anti-ID-Antigen-Bridging





## Overview of Design Elements

### 1. FDA recommends a multi-tiered approach

- **Tier 1:** A rapid, sensitive screening assay
- **Tier 2:** Positive should then be subjected to a confirmatory assay, such as ligand or antigen competition assay
- **Tier 3:** Further characterization: neutralization/class/isotype/titer


### 2. Aspects of Assay Development:

- Highly sensitive
- Able to detect all isotypes (IgM, IgE and IgG subtypes); carefully consider the avidity of control used to evaluate the assay; Should conduct assay performance test in the same concentration of matrix as that used to assess patient samples.


# Direct and Bridging ELISAs



## Direct Binding ELISA



## Bridging ELISA



✦ **Anti-ID (mAb or pAb) serve as the best positive controls** for IR assay development ✦

# Typical IR Assay Formats



## ◆ Direct and Bridging ELISA assays ★

- **Advantages:** Sensitive, inexpensive, equipment readily available
- **Disadvantages:** May not detect early immune response and may be influenced by high levels of circulating drug

## ◆ Radio-immuno precipitation assay (RIPA)

- **Advantages:** sensitive, inexpensive, equipment readily available
- **Disadvantages:** May not detect early immune response and may be influenced by high levels of circulating drug

## ◆ Surface Plasmon Resonance (SPR)

- **Advantages:** Method of choice for detecting early immune response and has Ab characterization capabilities
- **Disadvantages:** Expensive equipment, generally less sensitive than RIP or ELISA/ECL

## ◆ Electrochemiluminescence (ECL) assay

- **Advantages:** Sensitive, can be modified to respond in the presence of high levels of circulating drug
- **Disadvantages:** Equipment can be expensive, may not easily detect rapidly dissociating Abs


# Successful GenScript Case Studies



## Case study 1: Anti-ID mAb binds specifically to idiotype Ab, but not other negative controls

Clones	Antibody drug	Human IgG1	Normal IgG	Human serum
No.1	2.632	0.086	0.285	0.124
No.2	2.691	0.19	0.186	0.285
No.3	1.576	0.094	0.084	0.253
No.4	2.133	0.21	0.085	0.111
No.5	2.469	0.093	0.116	0.483
No.6	2.546	0.122	0.091	0.176
No.7	1.803	0.084	0.153	0.311
No.8	2.149	0.09	0.081	0.133
No.9	2.498	0.094	0.182	0.082
No.10	2.399	0.195	0.084	0.195

**Figure 1: Final sub-clones with high specificity to humanized Ab when screened against different counter screenings.**



**Figure 2: Final sub-clones with high specificity to humanized Ab when screened against isotype control, normal human IgG and human serum.**

# Successful GenScript Case Studies



**Case study 2 :** GenScript can generate high affinity anti-idiotypic antibody with sub-nanomolar level




Figure 3: Affinity Ranking of antibodies by using BIAcore; clones with sub-nM affinity

Clone 2




Figure 4: Affinity determination of Clone 2 antibodies by using BIAcore



## Case study 3 :





Figure 2: Specific inhibition of anti-Idiotypic mAbs in an antigen ligand blocking assay.



**GenScript has developed a panel of high affinity, high specificity anti-IDs for use in PK and IR assays**

## **Features of GenScript Anti-ID Platform:**

- High specificity and high affinity
- High speed production: **2-3 months**
- Antigen ligand blocking, epitope binning, and antibody pairing
- **Anti-ID polyclonal antibody packages** (rabbit & goat)
- **Anti-ID monoclonal antibody packages** (mouse & rat)
- Proven track record: **100% success rate**
- **PK** and **IR assay** development







# GenScript Anti-ID Packages



## Anti-Idiotypic Antibody Packages

Services	Anti-idiotypic mAb (SC1184)	Anti-idiotypic pAb (SC1185)
<b>Starting material</b>	Target antibody drug 2-3 mg	Target antibody drug 20 mg or more
<b>QC</b>	Cross-reactivity with control IgG <10%	Cross-reactivity with control IgG <10%
<b>Deliverables</b>	Hybridoma cell lines, supernatants and purified anti-ID antibody (optional)	0.5-3 mg purified anti-ID antibody/rabbit

## Ordering is easy:

Submit  
quotation  
request

A Technical  
Account Manager  
will e-mail you  
within 24 hours  
to finalize quote

Review and  
approve final  
order and price

To place order,  
provide credit  
card/PO  
information

Production of your  
project begins  
immediately

A Project Manager is in contact with  
you during the production process to  
give updates

<http://www.genscript.com/anti-idiotypic-antibody.html>

# Advantages of GenScript Custom Abs



- ◆ Sequence to purified antibody service with no need to provide an antigen

# Advantages of GenScript Custom Abs



- ◆ Sequence to purified antibody service with no need to provide an antigen
- ◆ Optimized immunization using our OptimumAntigen™ design tool and intelligent Antigen Strategy increasing specificity and affinity of antibodies

# Advantages of GenScript Custom Abs



- ◆ Sequence to purified antibody service with no need to provide an antigen
- ◆ Optimized immunization using our [OptimumAntigen™ design tool](#) and intelligent [Antigen Strategy](#) increasing specificity and affinity of antibodies
- ◆ Guaranteed results: quantity of antibodies or hybridoma, ELISA titer, and WB guarantee (varies with specific package)

# Advantages of GenScript Custom Abs



- ◆ Sequence to purified antibody service with no need to provide an antigen
- ◆ Optimized immunization using our [OptimumAntigen™ design tool](#) and intelligent [Antigen Strategy](#) increasing specificity and affinity of antibodies
- ◆ Guaranteed results: quantity of antibodies or hybridoma, ELISA titer, and WB guarantee (varies with specific package)
- ◆ Fast turnaround time: delivery of purified pAb or development of specific hybridoma in 45 days.
- ◆ Certified facility: AAALAC International accreditation and OLAW certification, demonstrating our commitment to responsible animal care and use.

# Variety of GenScript Antibody Services



## 🔬 Polyclonal Antibody Services

FAST pAb Services-PolyExpress™, Standard pAb Services

## 🔬 Phospho-Specific Antibody Services

Phospho-Specific pAb and mAb Services

## 🔬 Specialized Antibody Services

Antibody Drug Development, Immunoassay, Purification, Modifications, Conjugation

## 🔬 Monoclonal Antibody Services

FAST mAb Services-MonoExpress™, Custom mAb Services, Premium Hybridoma Services

## 🔬 Scale-up Antibody Services

Scale-up pAb and mAb Services, *In vivo* Ascite production, *In vitro* Roller bottle production







- Methylation protects microRNAs from an AGO1-associated activity that uridylates 5' RNA fragments generated by AGO1 cleavage.

*Yu B, Chen X, Vinovskis C, etc.*

*PNAS, (Apr 2014)*

- HYPERSENSITIVE TO HIGH LIGHT1 Interacts with LOW QUANTUM YIELD OF PHOTOSYSTEM II1 and Functions in Protection of Photosystem II from Photodamage in Arabidopsis.

*Wang HB, Wang J, Qi K, etc.*

*Plant Cell, (Mar 2014)*

- Tousled-like kinases phosphorylate Asf1 to promote histone supply during DNA replication.

*Groth A, Jensen ON, Nielsen ML, etc.*

*Nature Communications, (Mar 2014)*

- Dirigent domain-containing protein is part of the machinery required for formation of the lignin-based Casparian strip in the root.

*Hosmani PS, Kamiya T, Danku J, etc.*

*PNAS, (August 2013)*

- PfSETvs methylation of histone H3K36 represses virulence genes in Plasmodium falciparum



*Jiang L, Mu J, Zhang Q, Ni T, etc.*

*Nature, (July 2013)*

- Wheat Mds-1 encodes a heat-shock protein and governs susceptibility towards the Hessian fly gall midge.

*Liu X, Khajuria C, Li J, Trick HN, Huang L, etc*


*Nature Communications, (Jun 2013)*



# About GenScript



<http://www.genscript.com>








*Thank you for your participation  
We wish you all success in your research*

*Email me: [Michelle.Parker@GenScript.com](mailto:Michelle.Parker@GenScript.com)*

**Please complete the survey**

Register for other webinars in the GenScript Webinar Series @ <http://www.genscript.com/webinars.html>



November 25, 2014/8:00 am and 2:00 pm EST  
Peptide design strategy: basics, optimization, and  
application – *Tiffany Gupton Campolongo, PH.D.*



December 3, 2014/8:00 am and 2:00 pm EST  
Fusion partner for recombinant soluble protein  
production in *E. Coli* – *Keshav Vasanthavada*

# References



1. Karl Erik Hellstrom, Dale E. Yelton, H. Perry Fell, Donna Beaton, Margit Gayle, Michael Maclean, Maria Kahn, and Ingegerd Hellstrom. (1990) “**Epitope Mapping and use of Anti-Idiotypic Antibodies to the L6 Monoclonal Anticarcinoma**”. *Cancer Research* **50**, 2449-2454.
2. James A. Lofgren, Sripriya Dhandapani, Jason J. Pennucci, Christina M. Abbott, Daniel T. Mytych, Arunan Kaliyaperumal, Steven J. Swanson, and Michael C. Mullenix. (2007). “**Comparing ELISA and Surface Plasmon Resonance for Assessing Clinical Immunogenicity of Panitumumab**” *The Journal of Immunology* **178**: 7467-7472.
3. Karen Weiss, M.D., FDA. **News along the Pike**, December 3, 2003
4. Anthony R. Mire-Sluis, Yu Chen Barrett, Viswanath Devanarayan, Eugen Koren, Hank Liu, Mauricio Maia, Thomas Parish, George Scott, Gopi Shankar, Elizabeth Shores, Steven J. Swanson, Gary Taniguchi, Daniel Wierda, Linda A. Zuckerman (2004). “**Recommendations for the design and optimization of immunoassays used in the detection of host antibodies against biotechnology.**” *Journal of Immunological Methods* **289**: 1-16.
5. **Guidance for Industry Assay Development for Immunogenicity Testing of Therapeutic Protein** FDA, December 2009.
6. Jean W. Lee, Marian Kelley, Lindsay E. King, Jihong Yang, Hossein Salimi-Moosavi, Meina T. Tang, Jian-Feng Lu, John Kamerud, Ago Ahene, Heather Myler, Cindy Rogers. (2011) . “**Bioanalytical approaches to quantify “total” and “free” therapeutic antibodies and their targets: technical challenges and PK/PD applications over the course of drug development.**” *The AAPS Journal*, **13**(1): 99-110.
7. Shankar, G., Pendley, C., Stein, K. (2007). “**A risk-based bioanalytical strategy for the assessment of antibody immune responses against biological drugs.**” *Nature Biotechnology* **25**: 555-561.