

Avoiding Peptide Assay Failure: Hidden Problems and Solutions

Presented by:

Tiffany Gupton Campolongo, Ph.D.



Presentation Overview



- 1 Introduction
- 2 Key properties of custom peptides
- 3 Symptoms/Sequence evaluation
- 4 Solutions Summary & Resources
- 5 Q&A

About GenScript



Gene Services



Peptide Services



Protein Services



Antibody Services



Discovery Biology
Services



Catalog Products

Intrinsic Custom Peptide Properties



1

Contamination

- TFA counter-ions
- Protecting groups
- Truncated sequences
- Deletion sequences
 - Endotoxin
 - Water

2

Insolubility

- Precipitation
- Aggregation
- Secondary structure formation

3


Instability

- Oxidation
- Deliquescence
- Colorization
- Deamidation
- Hydrolysis
- Cyclization
- Secondary structure formation

Troubleshooting Peptide Assays




Evaluate Sequences



Examples:

- Calculate net charge
- Identify sequences having propensity for secondary structure formation
- Identify residues that bind contaminants

Identify Symptoms



Examples:


- Aberrant cellular or tissue growth
- Aberrant peptide bioactivity
- Abnormal physical peptide characteristics

Employ Solutions



- Pre-emptive sequence optimization
- Peptide preparation
- Assay optimization
- Specialized services


Overview of Solid Phase Peptide Synthesis



Quality Control Analyses




HPLC Trace Confirms peptide purity



Peak Table			
Peak#	Retention Time	Area	Area %
1	13.619	139296	1.453
2	13.916	359996	3.756
3	14.054	8707888	90.846
4	14.283	79644	0.831
5	14.400	37342	0.390
6	14.987	261173	2.725
Total		9585339	100.000

Contamination


MS Trace Confirms major ion species by mass





What is your custom peptide really made of?

Contamination




A peptide having 95% purity may actually only constitute: $(95\% \times 50\%) = 47.5\%$ of the lyophilized powder you receive

Note: Only an Amino Acid Test can tell you the exact concentration of your peptide

Trifluoroacetic acid (TFA) Counter-ion Contamination




Contamination



Symptoms

- Erratic cell or tissue viability, enzyme assay results
- IR spectroscopy data
- Reduced mass spec sensitivity
- Peptide degradation



Affected Sequences

All sequences, especially those containing positively charged:

- lysine
- histidine
- arginine




Solution

Protect your assays



TFA removal service

Case Study: TFA effects on cell and tissue culture



Effect of [TFA] on fetal rat osteoblast growth after 24 h.
 * $P < 0.05$ indicates data statistically significantly different from control.

Cell Type	Assay Type	[TFA] Effect
Fetal rat osteoblasts	Cell number	\leftrightarrow (10^{-9} M) \downarrow (10^{-7} M - 10^{-8} M)
Canine articular chondrocytes	Cell number	\leftrightarrow (4×10^{-9} M) \downarrow (4×10^{-8} M) \downarrow (4×10^{-7} M)
Neonatal mouse calvariae	Thymidine incorporation	\downarrow (2×10^{-7})

\uparrow Increase

\downarrow Decrease

\leftrightarrow No change

- TFA at concentrations higher than 10^{-9} M can decrease cell and tissue growth

*Figures adapted from Fig. 1 of Cornish J et al. Trifluoroacetate, a contaminant in purified proteins, inhibits proliferation of osteoblasts and chondrocytes. (1999) Am J Physiol Endocrinol Metab 277: E779-E783.

Case Study: TFA effects on cell and tissue culture





		Cell Type: Fetal Rat Osteoblasts					
		Amylin HCl-salt	Amylin TFA-salt	CGRP HCl-salt	CGRP TFA-salt	Rat calcitonin HCl-salt	Rat calcitonin TFA-salt
Parameter Measured	Cell number	↑	↔	↑↑	↑	↔	↓

		Tissue Type: Neonatal Mouse Calvariae	
		Amylin HCl-salt	Amylin TFA-salt
Parameter Measured	Thymidine incorporation	↑	↓

- Residual TFA in purified peptides can decrease cell and tissue growth

↑ Increase
 ↓ Decrease
 ↔ No change



Effect of TFA or HCl purified peptides on cell and tissue growth after 24 h


* $P < 0.05$ and ** $P < 0.003$ and indicates data statistically significantly different from control.

*Figures adapted from Fig. 2 and 5. of Cornish J et al. Trifluoroacetate, a contaminant in purified proteins, inhibits proliferation of osteoblasts and chondrocytes. (1999) Am J Physiol Endocrinol Metab 277: E779-E783.

GenScript's guaranteed TFA removal service



Contamination



Peptide Deliquescence




Contamination



Symptoms

- Peptide loses activity over time
- Lyophilized peptide converts from powder to crystal



Affected Sequences

Sequences containing:

- Ser
- Thr
- Lys
- Gly
- Arg



Solutions

Protect your peptide →

- Aliquot peptides
- Store peptides in a desiccator

*Deliquescence: Gradual dissolution by the absorption of moisture from the air

Case Study: Resolving deliquescence



Contamination




A - Lyophilized peptide after treatment with excipients

B - Lyophilized peptide without excipient treatment

Peptide Insolubility



Insolubility



Symptoms

Turbidity or precipitation
following peptide dissolution



Affected Sequences

Sequences:
>5 AA long, *and*
containing >50%:

- Tryptophan
- Isoleucine
- Leucine
- Phenylalanine
- Methionine
- Valine
- Tyrosine

Peptide insolubility solutions



Hydrophobicity

Peptide Protection

- Follow the solubility chart recommendations
- Request solubility testing



Solutions

Sequence Optimization

- Residue substitutions

Replace:

- Tryptophan
- Isoleucine
- Leucine
- Phenylalanine
- Methionine
- Valine
- Tyrosine



With:

- Aspartic acid
- Lysine
- Arginine
- Histidine
- Glutamic Acid


- Incorporate of O-acyl bonds
- Incorporate hydrophilic linkers (e.g. Lys-Lys-Lys-Lys)
- For libraries introduce a frame shift

Case study: O-acyl bond incorporation



Insolubility



β -amyloid solubility increased by O-acyl bond incorporation



	Click β -amyloid	Native β -amyloid
Solubility in water	15 mg/ml	0.14 mg/ml

Kiso et al. (2008) Click Peptides: Design and Applications

*GenScript's Click peptide service introduces a an o-acyl bond into the peptide sequence can greatly increase solubility



Example Solubility report results

Solvent	pH Value	Dissolution Achieved?	Peptide Concentration
Aqueous solution		Undissolved	N.A.
		Dissolved	mg/mL
		N.A.	N.A.
Buffer	pH 5	Dissolved	mg/mL
	pH 7	N.A.	N.A.
	pH 10	N.A.	N.A.
Organic solvent		Undissolved	N.A.
		Dissolved	mg/mL

Peptide Instability: Design Considerations




Instability

	Susceptible amino acids	Solutions
Cyclization	<ul style="list-style-type: none"> Asp-Gly, N-terminal Glu 	<ul style="list-style-type: none"> Substitute Asp with other amino acids Acetylate N-terminal Glu
Secondary structure formation	Multiple Glu, Ile, Leu, Phe, Thr, Tyr, Val	<ul style="list-style-type: none"> Asp for Glu Ser for Thr Pro or Gly every third residue
Oxidation	<ul style="list-style-type: none"> Cys Met 	<ul style="list-style-type: none"> Cysteine with Serine Methionine with Norleucine
Deliquescence	<ul style="list-style-type: none"> Ser Thr Lys Gly Arg 	Substitute non-essential charged residues with uncharged residues
Hydrolysis	<ul style="list-style-type: none"> Asp-Gly Asp-Pro Asp-Ser Asp-Gly 	Substitute Asp with other amino acids
Deletion sequences	Multiple adjacent Ser	Substitute Ser with other amino acids

Peptide Oxidation




Instability



Symptoms

- Peptide loses activity over time
- Peptide color changes from white to yellow, tan or brown



Affected Sequences

Sequences containing:

- Cysteine
- Tryptophan
- Methionine

Peptide Oxidation



Instability

Design Optimization



Solutions

Replace:

- Cysteine with Serine
- Methionine with Norleucine

Peptide Protection

- Aliquot peptides
- Flush peptides with argon gas and store in a tightly sealed container


Assay Protection

- Add reductants: DTT, TCEP, β -mercaptoethanol,
- Flush buffers with argon gas
- Reconstitute buffers and peptide with argon-flushed water in an anaerobic chamber
- Perform assays in anaerobic chamber

Case Study: Effects of Oxidation



Instability



*Argonshield™
Technology protects
peptides by
applying an argon
gas headspace and
airtight seal

Peptide after
synthesis and
purification

Peptide *without* ArgonShield™
protection after 6 months
of storage


(smaller peaks denote new
oxidized peptide species)

Peptides with
ArgonShield™ protection
after 6 months
of storage


Case study: Reducing/preserving oxidative residues



Instability



*Oxidation sensitive dipeptidyl transpeptidase enzyme activity assayed after experimental preparation and removal from an anaerobic chamber



*Oxidation sensitive dipeptidyl transpeptidase enzyme activity assayed with increasing levels of reductants

*Adapted from Fig. 4 and 5 of Gupton-Campolongo T et al. (2013) J. Phycol. 49:32-40

Peptide Color Change




Instability



Symptoms

- Erratic colorimetric assay data
- Lyophilized peptide is tan, yellow or brown



Affected Sequences

Sequences containing:

- Trp
- Tyr
- Cys
- Glu
- Asn



Solutions

Color removal service

*Lyophilized peptides can change color due to oxidation, or contamination by colored protective groups


Case Study: Removing Peptide Color



Instability


A

B




A

B



A

B



A - Lyophilized peptide after synthesis

B – Lyophilized peptide after proprietary de-colorization treatment

Resources: Stability Guidelines Chart



Do's and Dont's of Peptide Handling

For ALL peptides	<ul style="list-style-type: none">•DO aliquot lyophilized peptide according to daily experimental needs (Try GenScript's free aliquoting service)•DO avoid repeated freeze-thaw cycles•DON'T store peptides long-term in solution•DON'T repeatedly open the stock peptide vial
For peptides containing Cys, Met, or Trp residues	<ul style="list-style-type: none">•DO limit peptide exposure to air•DO purge assay buffers with argon or nitrogen gas•DO store peptides in tightly capped vials
For peptides containing Asp, Glu, Lys, Arg, or His	<ul style="list-style-type: none">•DO limit peptide exposure to air•DO store lyophilized peptides in a desiccator•DO store peptides in tightly capped vials
For peptides that must be stored in solution	<ul style="list-style-type: none">•DO avoid repeated freeze-thaw cycles•DO aliquot your peptide solution according to daily experimental needs•DO use sterile buffers to dissolve your peptide•DO filter your peptide using a 0.2 µm filter to remove bacterial contamination


http://www.genscript.com/peptide_storage_and_handling.html

Resources: Design Tools



Peptide MW calculator

- Calculates MW considering a selection of over **60** modifications



Peptide Property calculator

- Calculates isoelectric point
- Generates automatic solubility recommendations

- http://www.genscript.com/peptide_solubility_and_stability.html
- https://www.genscript.com/ssl-bin/peptide_mw
- https://www.genscript.com/ssl-bin/site2/peptide_calculation.cgi

Summary



Categories	Solutions	Specialized Services
Contamination	<ul style="list-style-type: none">▪ Amino acid substitutions▪ Storage in desiccator	<ul style="list-style-type: none">▪ TFA counter-ion removal▪ Moisture removal
Hydrophobicity	<ul style="list-style-type: none">▪ Amino acid substitutions▪ O-acyl bond incorporation▪ Hydrophilic linker incorporation▪ Frame shifting	<ul style="list-style-type: none">▪ O-acyl bond incorporation▪ Free Solubility testing
Instability	<ul style="list-style-type: none">▪ Amino acid substitutions▪ O-acyl bond incorporation▪ Aliquoting▪ Application of argon gas headspace▪ Use of reductants in assay▪ Perform assay in anaerobic chamber	<ul style="list-style-type: none">▪ Free Aliquoting▪ Argonsheid™▪ De-colorization service

GenScript's Specialized Services List



Service	Price	Turnaround Time (in addition to peptide synthesis)	Deliverables
Amino acid analysis	Typically \$100 + 5 mg additional peptide, sequence dependent	1 week	▪ Amino acid analysis test report
TFA Removal	Typically \$70 + 5 mg additional peptide, sequence dependent	2 weeks	▪ Lyophilized peptide ▪ TFA content test report
Click peptide service	Sequence dependent	No additional time	▪ Lyophilized peptide
Solubility testing	Free	No additional time	▪ Lyophilized peptide ▪ Solubility test report
ArgonShield™	Free	No additional time	▪ Lyophilized peptide in argon flushed, sealed vials
De-colorization treatment	Project dependent	1 week	▪ White Lyophilized peptide (by visual assessment)
Moisture removal	Project dependent	1 week	▪ Lyophilized peptide (by visual assessment)
Aliquoting	Free for first 5 vials, then \$2.00 per vial	No additional time	▪ Lyophilized peptide aliquoted into specified # of vials

Thank you for your participation
We wish you all success in your Research
Email me: Tiffany.Campolongo@GenScript.com



Register for other webinars in the GenScript Webinar Series @ <http://www.genscript.com/webinars.html>




May 21, 2014/ 11:00 am EST

Gene variant libraries: design, construction, and research applications - *Rachel Speer, Ph.D.*




May 28, 2014/ 1:00 pm EST

Protein or peptide antigen: choosing the optimal immunogen for antibody production - *Jessica Kaplunov, Ph.D.*




June 5, 2014/ 2:00 pm EST

Stem cell culture: choosing optimal conditions for expansion and differentiation - *Matthew Riolo, Ph.D.*



June 11, 2014/ 1:00 pm EST

Recombinant protein expression & purification: challenges and solutions - *Liyan Pang, Ph.D.*




June 18, 2014/ 2:00 pm EST

Can CRISPR/Cas9 off-target genomic editing be avoided? Ways to improve target specificity - *Maxine Chen, Ph.D.*



June 25, 2014/ 2:00 pm EST

Building a Synthetic Eukaryotic Genome – Sc2.0 - *Leslie Mitchell, Ph.D., NYU Langone Medical Center*



On Demand (Originally aired May 8, 2014)

Optimizing conditions for recombinant soluble protein production in *E. coli*, - *Keshav Vasanthavada*