

Version: 01

DATASHEET

Update: 08/09/2022

GenCRISPR™ LbuCas13a Nuclease

Cat. No.: Z03742

Product Introduction

GenCRISPR™ LbuCas13a Nuclease is an RNA-guided RNA endonuclease from *Leptotrichia buccalis*. Cas13a (previously known as C2c2) belongs to the Class 2 Type VI CRISPR-Cas system and contains a single protein effector with two HEPN domains. Cas13a exhibits two distinct ribonuclease activities. Cas13a possesses sequence-specific RNA cleavage activity *in vitro* and *in vivo*. And after specific cleavage, Cas13a activates collateral cleavage activity towards adjacent non-specific ssRNA sequences. The two distinct ribonuclease activities of Cas13a enable guide-RNA processing and RNA detection. The RNA detection empowered with Cas13a-crRNA complex and quenched fluorescent ssRNA as reporter has been reported to be versatile, rapid and high-specific.

Source: Recombinant Cas13a expressed by *E.coli*

Species: *Leptotrichia buccalis*

Accession#: C7NBY4

Tag: Tag-free

Apparent Molecular Weight: ~140 kDa, on SDS-PAGE under non-reducing conditions

Concentration: Please refer to the COA for the

specific lot.

Storage Buffer: 10 mM Tris-HCl, 300 mM NaCl, 1.0 mM DTT, 50% glycerol, pH 8.0.

Storage & Stability: Store at -20 °C for up to 12 months from the date of manufacture. Avoid repeated freeze-thaw cycles. **Do not store below -20 °C!**

Application: crRNA-dependent RNA cleavage

Quality Control Specifications

| Assay | Specifications |
|--------------------|---|
| Appearance | Clear, colorless liquid |
| Purity by SDS-PAGE | ≥ 90% |
| Activity | ≥ 85% by <i>in vitro</i> digestion of ssRNA |
| Residual DNase | Undetectable |
| Residual RNase | Undetectable |
| Endotoxin Level | < 0.2 EU/μg |

Reagents Supplied

| Components | Size | | | Storage |
|-------------------------------|--------|--------|--------|---------|
| GenCRISPR™ LbuCas13a Nuclease | 100 µg | 500 µg | 1 mg | -20 °C |
| 10 × Cas13a Reaction Buffer* | 1.5 ml | 1.5 ml | 1.5 ml | -20 °C |

* The reaction buffer is intended for Cas13a/crRNA-based Fluorescent Assays described in this manual, and not intended for other applications.

Protocol for Cas13a/crRNA-based Fluorescent Assays

1. Reagent preparation

- 1) Synthesize specific crRNA; we recommend GenScript's CRISPR Synthetic Guide RNA Services ([see www.genscript.com](http://www.genscript.com)) for crRNA synthesis.
- 2) Synthesize specific single-strand RNA (ssRNA) as substrate for cleavage; we recommend custom RNA synthesis using GenScript's RNA Synthesis Service ([see www.genscript.com](http://www.genscript.com)).
- 3) Synthesize quenched fluorescent ssRNA; we recommend custom RNA synthesis using GenScript RNA Synthesis Service ([see www.genscript.com](http://www.genscript.com)).
- 4) Prepare crRNA working solution at an appropriate concentration (e.g., 7.5 ng/µl) by diluting the stock solution with DEPC-Treated Water (nuclease-free) on ice.
- 5) Prepare specific ssRNA working solution at an appropriate concentration (e.g., 5 ng/µl) by diluting the stock solution with DEPC-Treated Water (nuclease-free) on ice.
- 6) The reaction volume is recommended to be set at 50 µl, and the reaction conditions should be optimized according to your specific applications.

2. Assay procedure

- 1) Assemble the reaction in a nuclease-free PCR strip tube for step 1 at room temperature on a clean bench as the following order in the table.
- 2) Assemble the reaction in a black ELISA Plate for step 2 at room temperature on a clean bench in the order indicated in the table.

| Step 1: Prepare the Cas13a/crRNA complex | |
|--|------------------|
| Components | Volume |
| 10 × Cas13a Reaction Buffer | 1 µl |
| 15 ng crRNA | 2 µl (7.5 ng/µl) |
| 200 ng GenCRISPR™ LbuCas13a Nuclease (Z03742) | 2 µl (100 ng/µl) |
| DEPC-Treated Water (nuclease-free) | 5 µl |
| Mix thoroughly, incubate for 10 min at 37 °C in a PCR instrument. | |
| Step 2: Fluorescent Assay | |
| Components | Volume |
| 10 ng ssRNA | 2 µl (5 ng/µl) |
| 10 × Cas13a Reaction Buffer | 4 µl |
| DEPC-Treated Water (nuclease-free) | 29 µl |
| Cas13a/crRNA complex (from step 1) | 10 µl |
| 10 pmol Fluorescent RNA reporter | 5 µl |
| Mix thoroughly, incubate for 30 min at 37 °C in an incubator. After incubation, read the | |

fluorescence intensity (Excitation wavelength: 494 nm; Emission wavelength: 518 nm) using a microplate reader.

References:

1. Abudayyeh, Omar O., et al. "C2c2 is a single-component programmable RNA-guided RNA-targeting CRISPR effector." *Science* 353.6299 (2016).
2. Gootenberg, Jonathan S., et al. "Nucleic acid detection with CRISPR-Cas13a/C2c2." *Science* 356.6336 (2017): 438-442.
3. Gootenberg, Jonathan S., et al. "Multiplexed and portable nucleic acid detection platform with Cas13, Cas12a, and Csm6." *Science* 360.6387 (2018): 439-444.
4. Iwasaki, Roman S., and Robert T. Batey. "SPRINT: a Cas13a-based platform for detection of small molecules." *Nucleic acids research* 48.17 (2020): e101-e101.
5. East-Seletsky, Alexandra, et al. "Two distinct RNase activities of CRISPR-C2c2 enable guide-RNA processing and RNA detection." *Nature* 538.7624 (2016): 270-273.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路 28 号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28 Yongxi Road, Jiangning District, Nanjing, Jiangsu, China