

# HEK293 Recombinant Promiscuous G Protein $G\alpha_{15}$ Cell Line



Technical Manual No. TM0299

Version 11242008

I	Introduction .....	1
II	Cell Line Information.....	1
III	Results.....	2
IV	References .....	3
V	Appendix .....	3
	User Agreement (Label License).....	6

## I. Introduction

The murine G-protein  $\alpha$ -subunit  $G\alpha_{15}$  belongs to the Gq protein subfamily, and it is also a promiscuous G protein because  $G\alpha_{15}$  can not only couple to Gq-coupled GPCRs but also to Gi-coupled GPCRs.  $G\alpha_{15}$  is widely used in GPCR-related research and pharmaceutical discoveries to change Gi-coupled GPCRs into Gq-coupled GPCRs, which are more easily detected by  $Ca^{2+}$  mobilization assays because when  $G\alpha_{15}$  is coupled by GPCRs, it can activate phospholipase C, increasing inositol trisphosphate, releasing  $Ca^{2+}$  from intracellular stores.

This manual describes establishment of a cell line and a protocol of pharmacologically validated for Gi-coupled GPCRs researches or Gi-coupled-GPCRs targeted high-throughput screening assay development. The cell line is stable through at least passage 15 and demonstrates a good dynamic range in fluorescent  $Ca^{2+}$  mobilization assays.

## II. Cell Line Information

- Catalog number: M00142
- Cell Line Name: HEK293/ $G\alpha_{15}$
- Description:


The  $G\alpha_{15}$  cDNA is amplified by PCR using a high fidelity enzyme and subcloned into the pIRES mammalian expression vector. The full-length ORF has been confirmed by sequencing. The  $G\alpha_{15}$  cell line is created by transfection of pIRES/ $G\alpha_{15}$  in HEK293 cells. The transfected cells are stably selected by 250  $\mu$ g/ml hygromycin. Single cell clones with high  $G\alpha_{15}$  inducibility and low background are isolated using ring cloning and serial dilution. The clones with largest dynamic range in  $Ca^{2+}$  mobilization are chosen for pharmacological and stability studies.
- Function: Cell based, functional assay for  $G\alpha_{15}$
- Quantity: 1 vial ( $2 \times 10^6$ ) frozen cells
- Passage Number Shipped: 2
- Stability Tested to: 16
- Host Cell: HEK293
- Cell Phenotype: Adherent/epithelial
- Antibiotics Selection: Hygromycin
- Recommended Storage: Liquid nitrogen, upon delivery
- Propagation Medium: DMEM, 10% FBS, 100  $\mu$ g/ml hygromycin, P/E
- Mycoplasma: Negative



### III. Results

Because  $G\alpha_{15}$  is only a G protein, we must transfect a Gi-coupled GPCR into HEK/ $G\alpha_{15}$  cells to verify that HEK293/ $G_{15}$  cells truly express functional  $G\alpha_{15}$  proteins.

#### Dopamine dose curve on HEK293/D2R/ $G\alpha_{15}$ cells



**Fig. 1.** Dose response of calcium currents as monitored with FlexStation plate reader upon treatment with dopamine. Assay was performed in accordance with the procedure described as below. Data represent means  $\pm$  SEM for triplicate samples (left one).  $EC_{50}$  value for dopamine dose response was determined using GraphPad Prism 4 software (right one).

**Note:** The HEK293/D2R/ $G\alpha_{15}$  stable cell line was established based on HEK293/ $G\alpha_{15}$  cells, which had been frozen and thawed twice.

### Assay Procedure

#### Fluorescent $Ca^{2+}$ -mobilization Assay of HEK293/ $G\alpha_{15}$ Cells

1. Lightly trypsinize cells and seed them into 96-well plate (Metrigel coating) at 25,000 cells/well (100  $\mu$ l per well) in complete culture medium.
2. Culture cells in 5%  $CO_2$  at 37°C for 24 hours (at least overnight).
3. Discard the growth media carefully using multi-channel pipette.
4. Wash cells once carefully with 100  $\mu$ l/well HBSS (containing 1 mM probenecid) using multi-channel pipette.
5. Stain cells with Fluo3/AM. Reconstitute the Fluo3/AM to a concentration of 2.0 mM for the stock solution. Dilute the Fluo3/AM of 2 mM stock to 4  $\mu$ M in HBSS as the loading solution (containing 1 mM probenecid, 1X red dye). Stain cells with 100  $\mu$ l/well of the loading solution.
6. Incubate the plate at 37°C in the dark for one hour.
7. Prepare agonist addition plates in advance of assay. To a 96-well plate, add 5X working concentration of Dopamine in HBSS, and add 25  $\mu$ l/well to cell plate.
8. Read with Flexstation using the specified settings and save data.

**Note:** All HBSS buffer used here containing 1 mM probenecid, including dye loading buffer, washing buffer and buffer to prepare agonist, etc. This solution is prepared fresh by diluting 200 mM probenecid stock solution (Dissolved in 200 mM NaOH) into HBSS Assay Buffer.



## IV. References

1. Gα<sub>15</sub> and Gα<sub>16</sub> couple a wide variety of receptors to phospholipase C. Offermanns S, and Simon M. JBC 270(25), 15175-15180, 1995.
2. Gz signaling: emerging divergence from Gi signaling. Ho M.K.C. and Wong Y.H. Oncogene 20, 1615-1625, 2001.

## V. Appendix

### Cell Culture Conditions

*Note: This clone often grows as clumps. It is necessary to separate cells into single cell suspensions whenever the cells are passaged or plated. MATRIGEL Matrix-coated flasks or plates are recommended for healthier cells.*

#### Complete Culture Medium:

DMEM: 90%, FBS: 10%, L-glutamine: 2.0 mM, Penicillinum: 100 µg/ml, Strep: 100 µg/ml, hygromycin: 150 µg/ml

#### Freezing Medium:

Growth medium plus 45% FBS and 10% DMSO

#### Thawing Cells:

1. Quickly thaw frozen cells in a 37°C water bath, agitating continuously.
2. Using a 1 ml pipette, slowly pipet the cells up and down five times and add, drop by drop, to a 15 ml centrifuge tube containing 5 ml of fresh prewarmed complete DMEM medium. Then centrifuge at 1,000 rpm for five minutes.
3. Discard the supernatant medium and resuspend the cell pellet in 5 ml of fresh prewarmed complete DMEM medium. Transfer cells to a T25 flask and incubate at 37°C with 5% CO<sub>2</sub> until the cells reach >90% confluence. The recovery rate for frozen cells is usually 90% or above.

#### Subculturing:

When the cells reach confluence, they need split. This cell line is normally split twice weekly at 1:8 to 1:15 dilutions.

1. Carefully aspirate all the media. Gently rinse the cell layer with appropriate amount of 0.2% trypsin-EDTA, and aspirate it off.
2. Wait for about 1-3 minutes. Then dislodge the cells by gently tapping the sides of flask or dish.
3. Resuspend cells with appropriate amount of complete DMEM medium, and split cells as desired.

#### Changing Medium:

This is normally done every other day.

1. Gently aspirate off medium.
2. Transfer fresh warm complete DMEM medium (37°C) into a flask (5 ml for T25 and 10 ml for T75).

#### Freezing Cells:

1. Repeat the steps 1-3 of subculturing section.
2. Centrifuge down the cells at 1,000 rpm for five minutes.



3. Aspirate off the supernatant and resuspend the cells in fresh freezing medium at a density of  $2-3 \times 10^6$  cells/ml. Add 1 ml cells per cryogenic vial.
4. Put the cryogenic vial of cells into cryo freezing container. Then transfer the container to a  $-80^\circ\text{C}$  environment and leave it there overnight.
5. Transfer cryogenic vial into liquid nitrogen ( $-196^\circ\text{C}$ ).

**Reagents & Consumables:**

1. DMEM: Dulbecco's Modified Eagle Medium powder, high-glucose (Gibco BRL, Cat #12100-046)
2. FBS: Fetal Bovine Serum (Hyclone, Cat #CH30160.03)
3. L-Glutamine: 200 mM (Gibco BRL, Cat # 25030-081)
4. Penicillin-Streptomycin, liquid (Gibco BRL, Cat # 15070-063)
5. Hygromycin B in PBS, 50 mg/ml (Invitrogen, Cat #10687-010)
6. 0.05% Trypsin-EDTA: (Gibco BRL, Cat # 25300)
7. DMSO: dimethyl sulphoxide, for molecular biology (Sigma, Cat #D8418)
8. Hepes: (Sigma Cat #H-3375)
9. Cryogenic Vial: (Corning Cat #430289)
10. Fluo-3 AM (Molecular Probes 28B2-4)
11. Probenecid: (Sigma Catalog #P 8761)
12. Dopamine: (Sigma Cat #H8502)
13. 96 Well Flat Clear Bottom White Polystyrene TC-Treated Microplates (Corning, #3610)

**Media and Solutions:**

1. PBS (for preparation of 500 ml)

1) KCl:	0.1 g
2) $\text{KH}_2\text{PO}_4$ :	0.1 g
3) NaCl:	4.0 g
4) $\text{Na}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$ :	1.4 g

Dissolve the above components in double-distilled water ( $\text{ddH}_2\text{O}$ ) and adjust pH to 7.4 with 0.1 N NaOH. Add  $\text{ddH}_2\text{O}$  to the final volume of 500 ml. Autoclave and store at  $4^\circ\text{C}$ .

2. Culture medium (for preparation of 1 L)


- 1) Measure out 950 ml distilled water to dissolve the media components, stirring gently until the solution becomes clear.
- 2) Add  $\text{NaHCO}_3$  3.7 g for high glucose DMEM.
- 3) Adjust pH of medium to 0.2-0.3 below the desired final working pH (using 1 N NaOH or 1 N HCL is recommended). Add slowly while stirring.
- 4) Dilute to 1 liter with  $\text{ddH}_2\text{O}$ .
- 5) Sterilize the medium immediately using the method of membrane filtration.  
Store at  $4^\circ\text{C}$ .

**Map of pIRES/  $\text{G}\alpha 15$** Name: pIRES/m $\text{G}\alpha 15$ Insert gene: Murine  $\text{G}\alpha 15$ 

Length of insert: 1128bp Sequence reference: Gene Bank M80632

Vector: pIRES (Clontech) Insert site : BamH I

Plasmid Map:

**Comments for pIRES/ $\text{G}\alpha 15$  6922 nucleotides**

CMV promoter: bases 309-810

 $\text{G}\alpha 15$ : bases 963-2090

Synthetic intron (IVS): bases 2136-2431

IRES: bases 2467-3069

Hygromycin resistance gene: bases 3067-4098

polyA signal: base 4391-4667

ColE1 origin of replication: base 5165-5784

Ampicillin resistance gene: bases 5917-6856

GenScript Corporation  
120 Centennial Ave., Piscataway, NJ 08854  
Tel: 732-885-9188, 732-885-9688  
Fax: 732-210-0262, 732-885-5878  
Email: [info@genscript.com](mailto:info@genscript.com)  
Web: <http://www.genscript.com>

**For Research Use Only.**



## **User Agreement (Label License) for HEK293 Recombinant Promiscuous G Protein Gα15 Cell Line**

**Catalog Number: M00142**

This is a legal agreement between you (Licensee) and GenScript Corporation governing use of GenScript's Stable GPCR cell line products and protocols provided to licensee. By using this product licensed to the licensee, the licensee is consenting to be legally bound by all of the terms of this user agreement. If you do not agree to all these terms, do not use the products or protocols. All products must be immediately returned to GenScript Corporation for a refund. The licensee will obtain no ownership interest in the product or any use/culture protocols accompanying the product other than through the perpetual limited license granted herein. If the product is licensed through an authorized GenScript distributor, the licensee will be obligated to disclose its identity to GenScript to ensure compliance with this user agreement.

All products sold by GenScript are for laboratory and animal research purposes only. The products are not to be used on humans, for consumption, or for other uses. The licensee may not sell or otherwise transfer this product or materials made using this product to a third party. The licensee may not use the GenScript products or materials for any commercial purposes. The vendee may transfer information or materials to a scientific collaborator only if such transfer is not for any commercial purpose. "Commercial purposes" means any activity by a vendee of this product for consideration and may include, but is not limited to the following: use of the product or its components in manufacturing; use of the product or its components for therapeutic, diagnostic, or prophylactic purposes; resale of the information or data developed from the product for research and development applications in a service business; or resale of the product or its components, whether or not such product or its components are resold for use in research.

GenScript retains ownership of all patents, copyrights, trademarks, trade secrets and other proprietary rights relating to or residing in the product and every portion thereof. The licensee agrees, (i) that it shall use all products solely in accordance with this agreement and that any such use of products will not violate any applicable law, regulation, judicial order, or injunction; and (ii) that it is not prohibited from receiving the products under U.S. export laws, that it is not a national of a country subject to U.S. trade sanctions, that it will not use the products in a location that is the subject of U.S. trade sanctions that would cover the products, and that, to its knowledge, it is not on the U.S. Department of Commerce's table of deny orders or is otherwise prohibited from obtaining goods of this sort from the United States.

For products that are subject to multiple limited use label licenses, the most restrictive terms apply. GenScript Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by GenScript Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine, or prophylactic product developed in research by the vendee in which this product or its components has been employed, provided that neither this product nor any of its components is used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Marketing Department, GenScript Corporation, 120 Centennial Avenue, Piscataway, New Jersey 08840. Phone: (732) 885-9188. Fax: (732) 210-0262. Email: [marketing@genscript.com](mailto:marketing@genscript.com)