


MS Spectrum



Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN

Sample Information Item
Date and Time : 2005-6-30 20:01:31
User : xie
Sample : YG-22-1
Inj. Volume : 1.000
MW : 2839.28

RP10424
PX00205509