


MS Spectrum



Sample Information Item : 2005-6-14 20:19:05
Date and Time : FJS
User : RP11277 AY-30-NH2
Sample : RP11277
Inj. Volume : 1.000
MW : 3200.62

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN