


MiniFresh[®] \ 'mi-nē 'fresh \ noun.

GenScript's fresh antibodies that are in small packages to help you reduce your research costs.



MiniFresh[®] Antibodies

MiniFresh[®] = \$50 antibodies

We provide small-package antibodies to keep your options open and your supplies fresh. Most of the antibodies are sold for only \$50 to reduce your research costs.

Small, because the MiniFresh[®] helps you to reduce your research costs;

Small, with big performance that helps you get quality results;

Small, because we concern about your costs and it is from GenScript.

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



MiniFresh® antibodies

Cat. No.	Name	Size	Price
A00131	Goat Anti-Rabbit IgG (H&L) Polyclonal Antibody	1 mg	\$38.50
A00132	Goat Anti-Chicken IgG (H&L) Polyclonal Antibody	1 mg	\$49.00
A00174	Rabbit Anti-His-tag Polyclonal Antibody	40 ug	\$50.00
A00170	Rabbit Anti-DYKDDDDK-tag Polyclonal Antibody	40 ug	\$50.00
A00172	Rabbit Anti-c-Myc-tag Polyclonal Antibody	40 ug	\$50.00
A00097	Rabbit Anti-GST-tag Polyclonal Antibody	40 ug	\$50.00
A00626	Rabbit Anti-Strep-tag II Polyclonal Antibody	40 ug	\$50.00
A00635	Rabbit Anti-CBP-tag Polyclonal Antibody	40 ug	\$50.00
A00191	Goat Anti-GAPDH Polyclonal Antibody	40 ug	\$50.00
A00683	Rabbit Anti-TAP-tag Polyclonal Antibody	40 ug	\$50.00
A00185	Mouse Anti-cGFP-tag Monoclonal Antibody	40 ug	\$50.00
A00180	Mouse Anti-Trx-tag Monoclonal Antibody	40 ug	\$50.00
A00682	Rabbit Anti-RFP-tag Polyclonal Antibody	40 ug	\$50.00
A00168	Goat Anti-HA-tag Polyclonal Antibody	40 ug	\$50.00
A00674	Rabbit Anti-Avi-tag Polyclonal Antibody	40 ug	\$50.00
A00623	Rabbit Anti-V5-tag Polyclonal Antibody	40 ug	\$50.00
A00190	Mouse Anti-MBP-tag Monoclonal Antibody	40 ug	\$50.00
A00636	Rabbit Anti-HAT-tag Polyclonal Antibody	40 ug	\$50.00
A00344	Rabbit Anti-EGFR (Ab-1070) Polyclonal Antibody	40 ug	\$50.00
A00044	Mouse Anti-Human IL-8 Monoclonal Antibody	100 ug	\$50.00
A00430	Rabbit Anti-HDAC5 (Ab-498) Polyclonal Antibody	40 ug	\$50.00
A00624	Rabbit Anti-HSV-tag Polyclonal Antibody	40 ug	\$50.00
A00447	Rabbit Anti-Cofilin (Ab-3) Polyclonal Antibody	40 ug	\$50.00
A00607	Rabbit Anti-MyoD (Ab-200) Polyclonal Antibody	40 ug	\$50.00
A00164	Mouse Anti-Thioredoxin (TRX) Monoclonal Antibody	40 ug	\$50.00
A00163	Mouse Anti-IP-10 Monoclonal Antibody	100 ug	\$50.00
A00316	Rabbit Anti-Estrogen Receptor- α (Ab-104) Polyclonal Antibody	40 ug	\$50.00
A00515	Rabbit Anti-CaMKII (Ab-286) Polyclonal Antibody	40 ug	\$50.00
A00380	Rabbit Anti-IGF-1R (Ab-1161) Polyclonal Antibody	40 ug	\$50.00
A00199	Rabbit Anti-VSV-G-tag Polyclonal Antibody	40 ug	\$50.00
A00300	Rabbit Anti-ATF-4 (Ab-245) Polyclonal Antibody	40 ug	\$50.00
A00518	Rabbit Anti-eIF2 α (Ab-51) Polyclonal Antibody	40 ug	\$50.00
A00517	Rabbit Anti-PKR (Ab-446) Polyclonal Antibody	40 ug	\$50.00
A00446	Rabbit Anti-cdc25A (Ab-75) Polyclonal Antibody	40 ug	\$50.00
A00684	Rabbit Anti-SNAP-tag Polyclonal Antibody	40 ug	\$50.00
A00279	Rabbit Anti-STAT4 (Ab-693) Polyclonal Antibody	40 ug	\$50.00
A00208	Rabbit Anti-PDK1 (Ab-241) Polyclonal Antibody	40 ug	\$50.00
A00888	Mouse Anti-Cre Recombinase Monoclonal Antibody	40 ug	\$50.00
A00580	Rabbit Anti-c-Jun (Ab-63) Polyclonal Antibody	40 ug	\$50.00
A00197	Rabbit Anti-c-Jun (Ab-73) Polyclonal Antibody	40 ug	\$50.00
A00508	Rabbit Anti-PKC θ (Ab-676) Polyclonal Antibody	40 ug	\$50.00
A00595	Rabbit Anti-p53 (Ab-15) Polyclonal Antibody	40 ug	\$50.00
A00129	Goat Anti-Trx-tag Polyclonal Antibody	40 ug	\$50.00
A00453	Rabbit Anti-I κ B- α (Ab-42) Polyclonal Antibody	40 ug	\$50.00
A00625	Rabbit Anti-S-tag Polyclonal Antibody	40 ug	\$50.00
A00586	Rabbit Anti-Pyk2 (Ab-402) Polyclonal Antibody	40 ug	\$50.00
A00610	Rabbit Anti-HSF1 (Ab-303) Polyclonal Antibody	40 ug	\$50.00
A00404	Rabbit Anti-G3BP-1 (Ab-232) Polyclonal Antibody	40 ug	\$50.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00583	Rabbit Anti-NF-κB p100/p52 (Ab-865) Polyclonal Antibody	40 ug	\$50.00
A00557	Rabbit Anti-SAPK/JNK (Ab-185) Polyclonal Antibody	40 ug	\$50.00
A00311	Rabbit Anti-BAD (Ab-155) Polyclonal Antibody	40 ug	\$50.00
A00669	Rabbit Anti-Ephrin-B2 (Ab-316) Polyclonal Antibody	40 ug	\$50.00
A00331	Rabbit Anti-HER2 (Ab-1248) Polyclonal Antibody	40 ug	\$50.00
A00352	Rabbit Anti-VEGFR2 (Ab-1175) Polyclonal Antibody	40 ug	\$50.00
A00553	Rabbit Anti-p38 MAPK (Ab-182) Polyclonal Antibody	40 ug	\$50.00
A00574	Rabbit Anti-eIF4E (Ab-209) Polyclonal Antibody	40 ug	\$50.00
A00413	Rabbit Anti-Caveolin-1 (Ab-14) Polyclonal Antibody	40 ug	\$50.00
A00573	Rabbit Anti-Stathmin 1 (Ab-15) Polyclonal Antibody	40 ug	\$50.00
A00303	Rabbit Anti-PTEN (Ab-370) Polyclonal Antibody	40 ug	\$50.00
A00325	Rabbit Anti-HER2 (Ab-877) Polyclonal Antibody	40 ug	\$50.00
A00302	Rabbit Anti-PTEN (Ab-380, 382, 383) Polyclonal Antibody	40 ug	\$50.00
A00388	Rabbit Anti-Chk2 (Ab-68) Polyclonal Antibody	40 ug	\$50.00
A00301	Rabbit Anti-Akt (Ab-473) Polyclonal Antibody	40 ug	\$50.00
A00189	Mouse Anti-HA-Tag Monoclonal Antibody	40 ug	\$50.00
A00590	Rabbit Anti-JunD (Ab-255) Polyclonal Antibody	40 ug	\$50.00
A00256	Rabbit Anti-JunB (Ab-259) Polyclonal Antibody	40 ug	\$50.00
A00429	Rabbit Anti-HDAC4 (Ab-632) Polyclonal Antibody	40 ug	\$50.00
A00450	Rabbit Anti-FKHRL1 (Ab-253) Polyclonal Antibody	40 ug	\$50.00
A00519	Rabbit Anti-CREB (Ab-129) Polyclonal Antibody	40 ug	\$50.00
A00565	Rabbit Anti-CDC2 (Ab-15) Polyclonal Antibody	40 ug	\$50.00
A00449	Rabbit Anti-Src (Ab-529) Polyclonal Antibody	40 ug	\$50.00
A00383	Rabbit Anti-IGF-1R (Ab-1165, 1166) Polyclonal Antibody	40 ug	\$50.00
A00725	Rabbit Anti-SRT-tag Polyclonal Antibody	40 ug	\$50.00
A00427	Rabbit Anti-FKHR (Ab-256) Polyclonal Antibody	40 ug	\$50.00
A00564	Rabbit Anti-p44/42 MAP Kinase (Ab-202) Polyclonal Antibody	40 ug	\$50.00
A00274	Rabbit Anti-STAT1 (Ab-701) Polyclonal Antibody	40 ug	\$50.00
A00619	Rabbit Anti-HRP (Horseradish Peroxidase) Polyclonal Antibody	200 ug	\$50.00
A00816	Rabbit Anti-PON2 Polyclonal Antibody	40 ug	\$50.00
A00424	Rabbit Anti-SEK1/MKK4 (Ab-80) Polyclonal Antibody	40 ug	\$50.00
A00672	Rabbit Anti-HSP90B (Ab-254) Polyclonal Antibody	40 ug	\$50.00
A00538	Rabbit Anti-EGFR (Ab-1110) Polyclonal Antibody	40 ug	\$50.00
A00701	Rabbit Anti-Serum Amyloid P Polyclonal Antibody	40 ug	\$50.00
A00228	Rabbit Anti-NF-κB p105/p50 (Ab-907) Polyclonal Antibody	40 ug	\$50.00
A00357	Rabbit Anti-VEGFR2 (Ab-951) Polyclonal Antibody	40 ug	\$50.00
A00582	Rabbit Anti-β-Catenin (Ab-37) Polyclonal Antibody	40 ug	\$50.00
A00511	Rabbit Anti-MARCKS (Ab-158) Polyclonal Antibody	40 ug	\$50.00
A00441	Rabbit Anti-Akt2 (Ab-474) Polyclonal Antibody	40 ug	\$50.00
A00310	Rabbit Anti-BAD (Ab-136) Polyclonal Antibody	40 ug	\$50.00
A00198	Rabbit Anti-Elk-1 (Ab-383) Polyclonal Antibody	40 ug	\$50.00
A00398	Rabbit Anti-Tau (Ab-205) Polyclonal Antibody	40 ug	\$50.00
A00642	Rabbit Anti-MEF2A (Ab-312) Polyclonal Antibody	40 ug	\$50.00
A00602	Rabbit Anti-CDK2 (Ab-160) Polyclonal Antibody	40 ug	\$50.00
A00196	Rabbit Anti-GSK3β (Ab-9) Polyclonal Antibody	40 ug	\$50.00
A00265	Rabbit Anti-Myc (Ab-373) Polyclonal Antibody	40 ug	\$50.00
A00222	Rabbit Anti-NF-κB p65 (Ab-536) Polyclonal Antibody	40 ug	\$50.00
A00307	Rabbit Anti-BCL-2 (Ab-70) Polyclonal Antibody	40 ug	\$50.00
A00529	Rabbit Anti-Histone H2A.X (Ab-139) Polyclonal Antibody	40 ug	\$50.00
A00550	Rabbit Anti-p70 S6 Kinase (Ab-421) Polyclonal Antibody	40 ug	\$50.00
A00412	Rabbit Anti-Rb (Ab-780) Polyclonal Antibody	40 ug	\$50.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00434	Rabbit Anti-FAK (Ab-925) Polyclonal Antibody	40 ug	\$50.00
A00411	Rabbit Anti-Rb (Ab-807) Polyclonal Antibody	40 ug	\$50.00
A00433	Rabbit Anti-cdc25C (Ab-216) Polyclonal Antibody	40 ug	\$50.00
A00658	Rabbit Anti-PLCy2 (Ab-753) Polyclonal Antibody	40 ug	\$50.00
A00409	Rabbit Anti-Paxillin (Ab-118) Polyclonal Antibody	40 ug	\$50.00
A00278	Rabbit Anti-STAT3 (Ab-727) Polyclonal Antibody	40 ug	\$50.00
A00451	Rabbit Anti-VASP (Ab-238) Polyclonal Antibody	40 ug	\$50.00
A00680	Rabbit Anti-Universal-tag Polyclonal Antibody	40 ug	\$50.00
A00655	Rabbit Anti-PKCβ (Ab-641) Polyclonal Antibody	40 ug	\$50.00
A00255	Rabbit Anti-JunB (Ab-79) Polyclonal Antibody	40 ug	\$50.00
A00588	Rabbit Anti-MSK1 (Ab-376) Polyclonal Antibody	40 ug	\$50.00
A00679	Rabbit Anti-S1-tag Polyclonal Antibody	40 ug	\$50.00
A00407	Rabbit Anti-β-Catenin (Ab-41, 45) Polyclonal Antibody	40 ug	\$50.00
A00587	Rabbit Anti-c-Jun (Ab-170) Polyclonal Antibody	40 ug	\$50.00
A00210	Rabbit Anti-GSK3α (Ab-21) Polyclonal Antibody	40 ug	\$50.00
A00209	Rabbit Anti-Raf-1 (Ab-259) Polyclonal Antibody	40 ug	\$50.00
A00706	Rabbit Anti-Prohibitin Polyclonal Antibody	40 ug	\$50.00
A00426	Rabbit Anti-Histone H3.1 (Ab-10) Polyclonal Antibody	40 ug	\$50.00
A00516	Rabbit Anti-p70 S6 Kinase (Ab-424) Polyclonal Antibody	40 ug	\$50.00
A00584	Rabbit Anti-Catenin (Ab-33) Polyclonal Antibody	40 ug	\$50.00
A00609	Rabbit Anti-IκB-α (Ab-32, 36) Polyclonal Antibody	40 ug	\$50.00
A00889	Mouse Anti-Enterokinase Monoclonal Antibody	40 ug	\$50.00
A00926	Rabbit Anti-FAK (Ab-397) Polyclonal Antibody	40 ug	\$50.00
A00513	Rabbit Anti-PKR (Ab-451) Polyclonal Antibody	40 ug	\$50.00
A00465	Rabbit Anti-S6 Ribosomal Protein (Ab-235) Polyclonal Antibody	40 ug	\$50.00
A00270	Rabbit Anti-GATA1 (Ab-142) Polyclonal Antibody	40 ug	\$50.00
A00670	Rabbit Anti-Ephrin-B2 (Ab-330) Polyclonal Antibody	40 ug	\$50.00
A00421	Rabbit Anti-IKKα (Ab-23) Polyclonal Antibody	40 ug	\$50.00
A00355	Rabbit Anti-VEGFR2 (Ab-1214) Polyclonal Antibody	40 ug	\$50.00
A00225	Rabbit Anti-NF-κB p100/p52 (Ab-869) Polyclonal Antibody	40 ug	\$50.00
A00533	Rabbit Anti-Merlin (Ab-518) Polyclonal Antibody	40 ug	\$50.00
A00603	Rabbit Anti-p90RSK (Ab-348) Polyclonal Antibody	40 ug	\$50.00
A00509	Rabbit Anti-PKCδ (Ab-645) Polyclonal Antibody	40 ug	\$50.00
A00578	Rabbit Anti-EGFR (Ab-1197) Polyclonal Antibody	40 ug	\$50.00
A00577	Rabbit Anti-EGFR (Ab-869) Polyclonal Antibody	40 ug	\$50.00
A00666	Rabbit Anti-EGFR (Ab-678) Polyclonal Antibody	40 ug	\$50.00
A00531	Rabbit Anti-Synapsin (Ab-9) Polyclonal Antibody	40 ug	\$50.00
A00286	Rabbit Anti-CREB (Ab-133) Polyclonal Antibody	40 ug	\$50.00
A00438	Rabbit Anti-MKK6 (Ab-207) Polyclonal Antibody	40 ug	\$50.00
A00459	Rabbit Anti-NF-κB p65 (Ab-529) Polyclonal Antibody	40 ug	\$50.00
A00264	Rabbit Anti-Myc (Ab-358) Polyclonal Antibody	40 ug	\$50.00
A00415	Rabbit Anti-MKK3 (Ab-189) Polyclonal Antibody	40 ug	\$50.00
A00350	Rabbit Anti-FAK (Ab-861) Polyclonal Antibody	40 ug	\$50.00
A00306	Rabbit Anti-BCL-2 (Ab-56) Polyclonal Antibody	40 ug	\$50.00
A00664	Rabbit Anti-AMPK1 (Ab-174) Polyclonal Antibody	40 ug	\$50.00
A00688	Rabbit Anti-β-Amyloid (1-40) Polyclonal Antibody	40 ug	\$50.00
A00436	Rabbit Anti-p27Kip1 (Ab-10) Polyclonal Antibody	40 ug	\$50.00
A00158	Goat Anti-Human IL-8 Polyclonal Antibody	40 ug	\$50.00
A00457	Rabbit Anti-Ezrin (Ab-566) Polyclonal Antibody	40 ug	\$50.00
A00327	Rabbit Anti-HER2 (Ab-1221,1222) Polyclonal Antibody	40 ug	\$50.00
A00918	Rabbit Anti-Beta-1,3-GalTase 2 Polyclonal Antibody	40 ug	\$50.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00304	Rabbit Anti-Akt (Ab-308) Polyclonal Antibody	40 ug	\$50.00
A00662	Rabbit Anti-14-3-3 ζ (Ab-58) Polyclonal Antibody	40 ug	\$50.00
A00456	Rabbit Anti-Paxillin (Ab-31) Polyclonal Antibody	40 ug	\$50.00
A00282	Rabbit Anti-STAT6 (Ab-641) Polyclonal Antibody	40 ug	\$50.00
A00661	Rabbit Anti-EGFR (Ab-1172) Polyclonal Antibody	40 ug	\$50.00
A00593	Rabbit Anti-Zap-70 (Ab-493) Polyclonal Antibody	40 ug	\$50.00
A00346	Rabbit Anti-EGFR (Ab-1092) Polyclonal Antibody	40 ug	\$50.00
A00949	Rabbit Anti-PKA RIIβ Polyclonal Antibody	40 ug	\$50.00
A00410	Rabbit Anti-Rb (Ab-795) Polyclonal Antibody	40 ug	\$50.00
A00547	Rabbit Anti-RelB (Ab-552) Polyclonal Antibody	40 ug	\$50.00
A00591	Rabbit Anti-DNA PKcs (Ab-2609) Polyclonal Antibody	40 ug	\$50.00
A00323	Rabbit Anti-Estrogen Receptor-α (Ab-167) Polyclonal Antibody	40 ug	\$50.00
A00628	Rabbit Anti-AU1-tag Polyclonal Antibody	40 ug	\$50.00
A00521	Rabbit Anti-Cortactin (Ab-421) Polyclonal Antibody	40 ug	\$50.00
A00386	Rabbit Anti-Integrin β3 (Ab-773) Polyclonal Antibody	40 ug	\$50.00
A00213	Rabbit Anti-NF-κB p65 (Ab-435) Polyclonal Antibody	40 ug	\$50.00
A00545	Rabbit Anti-Tau (Ab-212) Polyclonal Antibody	40 ug	\$50.00
A00235	Rabbit Anti-c-Jun (Ab-239) Polyclonal Antibody	40 ug	\$50.00
A00656	Rabbit Anti-PKCθ (Ab-695) Polyclonal Antibody	40 ug	\$50.00
A00184	Mouse Anti-M-CSF Monoclonal Antibody	100 ug	\$50.00
A00589	Rabbit Anti-Bcr (Ab-177) Polyclonal Antibody	40 ug	\$50.00
A00212	Rabbit Anti-NF-κB p65 (Ab-276) Polyclonal Antibody	40 ug	\$50.00
A00931	Rabbit Anti-α-Catenin (Ab-641) Polyclonal Antibody	40 ug	\$50.00
A00211	Rabbit Anti-MEK2 (Ab-394) Polyclonal Antibody	40 ug	\$50.00
A01046	Rabbit Anti-GAP43 (Ab-41) Polyclonal Antibody	40 ug	\$50.00
A00276	Rabbit Anti-STAT3 (Ab-705) Polyclonal Antibody	40 ug	\$50.00
A00428	Rabbit Anti-HDAC2 (Ab-394) Polyclonal Antibody	40 ug	\$50.00
A00903	Rabbit Anti-PKC β II Polyclonal Antibody	40 ug	\$50.00
A00233	Rabbit Anti-c-Jun (Ab-93) Polyclonal Antibody	40 ug	\$50.00
A00930	Rabbit Anti-Shc1 (Ab-427) Polyclonal Antibody	40 ug	\$50.00
A00677	Rabbit Anti-tag•100-tag Polyclonal Antibody	40 ug	\$50.00
A01048	Rabbit Anti-Tau (Ser356) Polyclonal Antibody	40 ug	\$50.00
A00542	Rabbit Anti-NF-κB p65 (Ab-311) Polyclonal Antibody	40 ug	\$50.00
A00893	Rabbit Anti-AMPK Alpha-1 Chain Polyclonal Antibody	40 ug	\$50.00
A00448	Rabbit Anti-Vav (Ab-174) Polyclonal Antibody	40 ug	\$50.00
A00929	Rabbit Anti-NFκB-p105/p50 (Ab-927) Polyclonal Antibody	40 ug	\$50.00
A00317	Rabbit Anti-Estrogen Receptor-α (Ab-106) Polyclonal Antibody	40 ug	\$50.00
A00676	Rabbit Anti-DDDDK-tag Polyclonal Antibody	40 ug	\$50.00
A01049	Rabbit Anti-Tau (Ser400) Polyclonal Antibody	40 ug	\$50.00
A00892	Rabbit Anti-Connexin 30 Polyclonal Antibody	40 ug	\$50.00
A00231	Rabbit Anti-c-Jun (Ab-91) Polyclonal Antibody	40 ug	\$50.00
A00651	Rabbit Anti-PAK1 (Ab-212) Polyclonal Antibody	40 ug	\$50.00
A00928	Rabbit Anti-IκB-β (Ab-23) Polyclonal Antibody	40 ug	\$50.00
A00585	Rabbit Anti-NF-κB p105/p50 (Ab-337) Polyclonal Antibody	40 ug	\$50.00
A00675	Rabbit Anti-PKC ε Polyclonal Antibody	40 ug	\$50.00
A01050	Rabbit Anti-Tau (Ser409) Polyclonal Antibody	40 ug	\$50.00
A00273	Rabbit Anti-GATA1 (Ab-310) Polyclonal Antibody	40 ug	\$50.00
A00425	Rabbit Anti-ASK1 (Ab-966) Polyclonal Antibody	40 ug	\$50.00
A00891	Rabbit Anti-Cytochrome P450 2C23 Polyclonal Antibody	40 ug	\$50.00
A00927	Rabbit Anti-IGF-1R (Ab-1346) Polyclonal Antibody	40 ug	\$50.00
A01012	Mouse Anti-WISP2 Monoclonal Antibody	40 ul	\$50.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00560	Rabbit Anti-HSP27 (Ab-78) Polyclonal Antibody	40 ug	\$50.00
A00445	Rabbit Anti-AFX (Ab-197) Polyclonal Antibody	40 ug	\$50.00
A00514	Rabbit Anti-BIM (Ab-65/69) Polyclonal Antibody	40 ug	\$50.00
A00973	Rabbit Anti-Runx3 Polyclonal Antibody	40 ug	\$50.00
A00608	Rabbit Anti-Chk1 (Ab-317) Polyclonal Antibody	40 ug	\$50.00
A00423	Rabbit Anti-SEK1/MKK4 (Ab-261) Polyclonal Antibody	40 ug	\$50.00
A00444	Rabbit Anti-FKHR (Ab-319) Polyclonal Antibody	40 ug	\$50.00
A00925	Rabbit Anti-ADD1 (Ab-726) Polyclonal Antibody	40 ug	\$50.00
A00671	Rabbit Anti-HSF1 (Ab-307) Polyclonal Antibody	40 ug	\$50.00
A00972	Rabbit Anti-RAGE Polyclonal Antibody	40 ug	\$50.00
A00700	Rabbit Anti-APP FM (339-348) Polyclonal Antibody	40 ug	\$50.00
A00401	Rabbit Anti-Tau (Ab-404) Polyclonal Antibody	40 ug	\$50.00
A00422	Rabbit Anti-ASK1 (Ab-83) Polyclonal Antibody	40 ug	\$50.00
A00887	Rabbit Anti-IDO Polyclonal Antibody	40 ug	\$50.00
A00227	Rabbit Anti-NF-κB p105/p50 (Ab-893) Polyclonal Antibody	40 ug	\$50.00
A00558	Rabbit Anti-HSP27 (Ab-82) Polyclonal Antibody	40 ug	\$50.00
A00512	Rabbit Anti-Gab1 (Ab-627) Polyclonal Antibody	40 ug	\$50.00
A00581	Rabbit Anti-NF-κB p65 (Ab-254) Polyclonal Antibody	40 ug	\$50.00
A00464	Rabbit Anti-Met (Ab-1234) Polyclonal Antibody	40 ug	\$50.00
A00970	Rabbit Anti-SHP-1 (Ab-536) Polyclonal Antibody	40 ug	\$50.00
A00606	Rabbit Anti-PKD/PKCμ (Ab-738) Polyclonal Antibody	40 ug	\$50.00
A00400	Rabbit Anti-Tau (Ab-231) Polyclonal Antibody	40 ug	\$50.00
A00442	Rabbit Anti-c-Abl (Ab-412) Polyclonal Antibody	40 ug	\$50.00
A00463	Rabbit Anti-Stathmin 1 (Ab-37) Polyclonal Antibody	40 ug	\$50.00
A00268	Rabbit Anti-MEF2A (Ab-319) Polyclonal Antibody	40 ug	\$50.00
A00605	Rabbit Anti-PKD/PKCμ (Ab-910) Polyclonal Antibody	40 ug	\$50.00
A00696	Rabbit Anti-Amyloid Precursor Protein (747-761) Polyclonal Antibody	40 ug	\$50.00
A00885	Rabbit Anti-α-Actin-1 Polyclonal Antibody	40 ug	\$50.00
A00510	Rabbit Anti-MEK1 (Ab-291) Polyclonal Antibody	40 ug	\$50.00
A00462	Rabbit Anti-Stathmin 1 (Ab-24) Polyclonal Antibody	40 ug	\$50.00
A00971	Rabbit Anti-LC3A/Apg-8 Polyclonal Antibody	40 ug	\$50.00
A00695	Rabbit Anti-Amyloid Precursor Protein (671-685) Polyclonal Antibody	40 ug	\$50.00
A00419	Rabbit Anti-JAK2 (Ab-221) Polyclonal Antibody	40 ug	\$50.00
A00884	Rabbit Anti-APH1 Polyclonal Antibody	40 ug	\$50.00
A00555	Rabbit Anti-NF-κB p105/p50 (Ab-932) Polyclonal Antibody	40 ug	\$50.00
A00309	Rabbit Anti-BAD (Ab-112) Polyclonal Antibody	40 ug	\$50.00
A00461	Rabbit Anti-4E-BP1 (Ab-45) Polyclonal Antibody	40 ug	\$50.00
A00966	Rabbit Anti-Amyloid Precursor Protein (Ab-730) Polyclonal Antibody	40 ug	\$50.00
A00266	Rabbit Anti-Elk-1 (Ab-389) Polyclonal Antibody	40 ug	\$50.00
A00694	Rabbit Anti-Amyloid Precursor Protein (653-662) Polyclonal Antibody	40 ug	\$50.00
A00397	Rabbit Anti-Tau (Ab-181) Polyclonal Antibody	40 ug	\$50.00
A00418	Rabbit Anti-JAK1 (Ab-1022) Polyclonal Antibody	40 ug	\$50.00
A00883	Rabbit Anti-Nicastrin Polyclonal Antibody	40 ug	\$50.00
A00638	Rabbit Anti-B-tag Polyclonal Antibody	40 ug	\$50.00
A00308	Rabbit Anti-BCL-XL (Ab-62) Polyclonal Antibody	40 ug	\$50.00
A00460	Rabbit Anti-4E-BP1 (Ab-36) Polyclonal Antibody	40 ug	\$50.00
A00601	Rabbit Anti-p90RSK (Ab-352) Polyclonal Antibody	40 ug	\$50.00
A00693	Rabbit Anti-Amyloid Precursor Protein (292-312) Polyclonal Antibody	40 ug	\$50.00
A00416	Rabbit Anti-TYK2 (Ab-1054) Polyclonal Antibody	40 ug	\$50.00
A00881	Rabbit Anti-Presenilin 1 Polyclonal Antibody	40 ug	\$50.00
A00637	Rabbit Anti-Protein C (HPC4)-tag Polyclonal Antibody	40 ug	\$50.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00963	Rabbit Anti-Piwi-like Protein 1 Polyclonal Antibody	40 ug	\$50.00
A00600	Rabbit Anti-Tau (Ab-262) Polyclonal Antibody	40 ug	\$50.00
A00692	Rabbit Anti-Amyloid Precursor Protein (44-63) Polyclonal Antibody	40 ug	\$50.00
A00880	Rabbit Anti-BACE1 (485-501) Polyclonal Antibody	40 ug	\$50.00
A00938	Rabbit Anti-V-ATPase Subunit A Polyclonal Antibody	40 ug	\$50.00
A00575	Rabbit Anti-IRS-1 (Ab-636) Polyclonal Antibody	40 ug	\$50.00
A00458	Rabbit Anti-VASP (Ab-157) Polyclonal Antibody	40 ug	\$50.00
A00962	Rabbit Anti-Peroxiredoxin-6 Polyclonal Antibody	40 ug	\$50.00
A00598	Rabbit Anti-p53 (Ab-315) Polyclonal Antibody	40 ug	\$50.00
A00414	Rabbit Anti-Src (Ab-418) Polyclonal Antibody	40 ug	\$50.00
A00879	Rabbit Anti-BACE1 (46-65) Polyclonal Antibody	40 ug	\$50.00
A01038	Rabbit Anti-ICAM-1 (Ab-512) Polyclonal Antibody	40 ug	\$50.00
A00284	Rabbit Anti-STAT6 (Ab-645) Polyclonal Antibody	40 ug	\$50.00
A01005	Rabbit Anti-Superoxide Dismutase 1 Polyclonal Antibody	40 ug	\$50.00
A00305	Rabbit Anti-p95/NBS1 (Ab-343) Polyclonal Antibody	40 ug	\$50.00
A00961	Rabbit Anti-PGK1 Polyclonal Antibody	40 ug	\$50.00
A00262	Rabbit Anti-Myc (Ab-58) Polyclonal Antibody	40 ug	\$50.00
A00596	Rabbit Anti-p53 (Ab-46) Polyclonal Antibody	40 ug	\$50.00
A00687	Rabbit Anti-β-Amyloid (25-35) Polyclonal Antibody	40 ug	\$50.00
A00393	Rabbit Anti-Tau (Ab-235) Polyclonal Antibody	40 ug	\$50.00
A00878	Rabbit Anti-BACE1 (26-45) Polyclonal Antibody	40 ug	\$50.00
A00219	Rabbit Anti-NF-κB p65 (Ab-468) Polyclonal Antibody	40 ug	\$50.00
A01039	Rabbit Anti-PLCγ1 (Ab-783) Polyclonal Antibody	40 ug	\$50.00
A00901	Rabbit Anti-PKCα Polyclonal Antibody	40 ug	\$50.00
A00348	Rabbit Anti-p73 (Ab-99) Polyclonal Antibody	40 ug	\$50.00
A00960	Rabbit Anti-V-ATPase Subunit B 2 Polyclonal Antibody	40 ug	\$50.00
A00261	Rabbit Anti-ATF-2 (Ab-94/112) Polyclonal Antibody	40 ug	\$50.00
A00685	Rabbit Anti-β-Amyloid (1-17) Polyclonal Antibody	40 ug	\$50.00
A00392	Rabbit Anti-Ezrin (Ab-353) Polyclonal Antibody	40 ug	\$50.00
A01040	Rabbit Anti-PAK1/PAK2/PAK3 (Ab-423/402/421) Polyclonal Antibody	40 ug	\$50.00
A00632	Rabbit Anti-KT3-tag Polyclonal Antibody	40 ug	\$50.00
A00939	Rabbit Anti-Estrogen Receptor Beta Polyclonal Antibody	40 ug	\$50.00
A00455	Rabbit Anti-SMC1 (Ab-957) Polyclonal Antibody	40 ug	\$50.00
A00954	Rabbit Anti-ApoM Polyclonal Antibody	40 ug	\$50.00
A00260	Rabbit Anti-ATF-2 (Ab-55/73) Polyclonal Antibody	40 ug	\$50.00
A00389	Rabbit Anti-Tau (Ab-396) Polyclonal Antibody	40 ug	\$50.00
A01041	Rabbit Anti-MEK1 (Ab-217) Polyclonal Antibody	40 ug	\$50.00
A00281	Rabbit Anti-STAT5A (Ab-780) Polyclonal Antibody	40 ug	\$50.00
A00631	Rabbit Anti-E-tag Polyclonal Antibody	40 ug	\$50.00
A00916	Rabbit Anti-Mitoferrin-1 Polyclonal Antibody	40 ug	\$50.00
A00238	Rabbit Anti-c-Jun (Ab-243) Polyclonal Antibody	40 ug	\$50.00
A00660	Rabbit Anti-Rac1 (Ab-71) Polyclonal Antibody	40 ug	\$50.00
A00454	Rabbit Anti-STAT1 (Ab-727) Polyclonal Antibody	40 ug	\$50.00
A00259	Rabbit Anti-ATF-2 (Ab-53/71) Polyclonal Antibody	40 ug	\$50.00
A00592	Rabbit Anti-Elk-1 (Ab-417) Polyclonal Antibody	40 ug	\$50.00
A01042	Rabbit Anti-IRS-1 (Ab-639) Polyclonal Antibody	40 ug	\$50.00
A00280	Rabbit Anti-STAT5A (Ab-694) Polyclonal Antibody	40 ug	\$50.00
A00630	Rabbit Anti-E2-tag Polyclonal Antibody	40 ug	\$50.00
A00432	Rabbit Anti-HSP27 (Ab-15) Polyclonal Antibody	40 ug	\$50.00
A00569	Rabbit Anti-Met (Ab-1349) Polyclonal Antibody	40 ug	\$50.00
A00942	Rabbit Anti-Tau (Ab-217) Polyclonal Antibody	40 ug	\$50.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00258	Rabbit Anti-ATF-2 (Ab-51/69) Polyclonal Antibody	40 ug	\$50.00
A01043	Rabbit Anti-GluR1 (Ab-849) Polyclonal Antibody	40 ug	\$50.00
A00431	Rabbit Anti-HDAC8 (Ab-39) Polyclonal Antibody	40 ug	\$50.00
A00452	Rabbit Anti-MEK1 (Ab-221) Polyclonal Antibody	40 ug	\$50.00
A00933	Rabbit Anti-Niban Polyclonal Antibody	40 ug	\$50.00
A00257	Rabbit Anti-ATF-2 (Ab-44/62) Polyclonal Antibody	40 ug	\$50.00
A00322	Rabbit Anti-Estrogen Receptor-α (Ab-118) Polyclonal Antibody	40 ug	\$50.00
A00681	Rabbit Anti-MAT-tag Polyclonal Antibody	40 ug	\$50.00
A01044	Rabbit Anti-p62Dok (Ab-362) Polyclonal Antibody	40 ug	\$50.00
A00908	Rabbit Anti-PKC ζ Polyclonal Antibody	40 ug	\$50.00
A00567	Rabbit Anti-BRCA1 (Ab-1423) Polyclonal Antibody	40 ug	\$50.00
A00932	Rabbit Anti-CARM1 (Ab-228) Polyclonal Antibody	40 ug	\$50.00
A00520	Rabbit Anti-Cortactin (Ab-466) Polyclonal Antibody	40 ug	\$50.00
A01045	Rabbit Anti-p62Dok (Ab-398) Polyclonal Antibody	40 ug	\$50.00
A00544	Rabbit Anti-Connexin 43 (Ab-367) Polyclonal Antibody	40 ug	\$50.00
A00905	Rabbit Anti-PKCδ Polyclonal Antibody	40 ug	\$50.00
A00852	Mouse Anti-OVA Monoclonal Antibody	40 ug	\$55.00
A00203	Goat Anti-HA-tag [Biotin] Polyclonal Antibody	40 ug	\$55.00
A00130	Rabbit Anti-GST-tag [HRP] Polyclonal Antibody	40 ug	\$55.00
A00169	Goat Anti-HA-tag [HRP] Polyclonal Antibody	40 ug	\$55.00
A00173	Rabbit Anti-c-Myc-tag [HRP] Polyclonal Antibody	40 ug	\$55.00
A00886	Mouse Anti-GFP [HRP] Monoclonal Antibody	40 ug	\$55.00
A00877	Rabbit Anti-V5-tag [HRP] Polyclonal Antibody	40 ug	\$55.00
A00202	Rabbit Anti-GST-tag [Biotin] Polyclonal Antibody	40 ug	\$55.00
A00846	Mouse Anti-LAL Monoclonal Antibody	40 ug	\$55.00
A00837	Mouse Anti-BSA Monoclonal Antibody	40 ug	\$55.00
A00815	Mouse Anti-CIB1 Monoclonal Antibody	40 ug	\$55.00
A00847	Mouse Anti-GLP Monoclonal Antibody	40 ug	\$55.00
A00737	Mouse Anti-Visfatin Monoclonal Antibody	40 ug	\$55.00
A00818	Mouse Anti-MAPK10 Monoclonal Antibody	40 ug	\$55.00
A00808	Mouse Anti-Dynamin-1 Monoclonal Antibody	40 ug	\$55.00
A00734	Mouse Anti-MCP-1 Monoclonal Antibody	40 ug	\$55.00
A00824	Mouse Anti-Vimentin Monoclonal Antibody	40 ug	\$55.00
A00872	Rabbit Anti-VSV-G-tag [HRP] Polyclonal Antibody	40 ug	\$55.00
A00819	Mouse Anti-Desmin Monoclonal Antibody	40 ug	\$55.00
A00817	Mouse Anti-GSTP1 Monoclonal Antibody	40 ug	\$55.00
A00814	Mouse Anti-CK1 Monoclonal Antibody	40 ug	\$55.00
A00204	Rabbit Anti-c-Myc-tag [Biotin] Polyclonal Antibody	40 ug	\$55.00
A00850	Mouse Anti-ER Monoclonal Antibody	40 ug	\$55.00
A00812	Mouse Anti-SNCG Monoclonal Antibody	40 ug	\$55.00
A00811	Mouse Anti-SNCA Monoclonal Antibody	40 ug	\$55.00
A00853	Mouse Anti-TUG Monoclonal Antibody	40 ug	\$55.00
A00839	Mouse Anti-GAPDH Monoclonal Antibody	40 ug	\$55.00
A00749	Mouse Anti-EphB6 Monoclonal Antibody	40 ug	\$55.00
A00838	Mouse Anti-Clenbuterol Monoclonal Antibody	40 ug	\$55.00
A00806	Mouse Anti-APP Monoclonal Antibody	40 ug	\$55.00
A00742	Mouse Anti-Human CD19 Monoclonal Antibody	40 ug	\$55.00
A00741	Mouse Anti-Human CD14 Monoclonal Antibody	40 ug	\$55.00
A00799	Mouse Anti-MSH2 Monoclonal Antibody	40 ug	\$55.00
A00796	Mouse Anti-Rab25 Monoclonal Antibody	40 ug	\$55.00
A00829	Mouse Anti-CD31 Monoclonal Antibody	40 ug	\$55.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00828	Mouse Anti-P16 (Mouse and Human) Monoclonal Antibody	40 ug	\$55.00
A00731	Mouse Anti-IFN-γ Monoclonal Antibody	40 ug	\$55.00
A00787	Mouse Anti-TIP60 Monoclonal Antibody	40 ug	\$55.00
A00783	Mouse Anti-BTK Monoclonal Antibody	40 ug	\$55.00
A00782	Mouse Anti-Lyn Monoclonal Antibody	40 ug	\$55.00
A00781	Mouse Anti-PTK6 Monoclonal Antibody	40 ug	\$55.00
A00780	Mouse Anti-AURKB Monoclonal Antibody	40 ug	\$55.00
A00778	Mouse Anti-ICKB Monoclonal Antibody	40 ug	\$55.00
A00777	Mouse Anti-GRN Monoclonal Antibody	40 ug	\$55.00
A00813	Mouse Anti-Lyn Monoclonal Antibody	40 ug	\$55.00
A00776	Mouse Anti-Dynamin-2 Monoclonal Antibody	40 ug	\$55.00
A00763	Mouse Anti-Tyro3 Monoclonal Antibody	40 ug	\$55.00
A00809	Mouse Anti-SRA Monoclonal Antibody	40 ug	\$55.00
A00744	Mouse Anti-Foxp3 Monoclonal Antibody	40 ug	\$55.00
A00804	Mouse Anti-Calcyclin Monoclonal Antibody	40 ug	\$55.00
A00835	Mouse Anti-Apoa5 Monoclonal Antibody	40 ug	\$55.00
A00802	Mouse Anti-LSD1 Monoclonal Antibody	40 ug	\$55.00
A00832	Mouse Anti-IGF1R Monoclonal Antibody	40 ug	\$55.00
A00831	Mouse Anti-Human P16 Monoclonal Antibody	40 ug	\$55.00
A00740	Mouse Anti-Human CD10 Monoclonal Antibody	40 ug	\$55.00
A00830	Mouse Anti-RSK1 Monoclonal Antibody	40 ug	\$55.00
A00794	Mouse Anti-Bcl-10 Monoclonal Antibody	40 ug	\$55.00
A00735	Mouse Anti-VEGI polyclonal antibody	40 ug	\$55.00
A00793	Mouse Anti-FES Monoclonal Antibody	40 ug	\$55.00
A00876	Rabbit Anti-E2-tag [HRP] polyclonal antibody	40 ug	\$55.00
A00792	Mouse Anti-PRAK Monoclonal Antibody	40 ug	\$55.00
A00827	Mouse Anti-P16 (Mouse/Human) Monoclonal Antibody	40 ug	\$55.00
A00874	Rabbit Anti-E-tag [HRP] Polyclonal Antibody	40 ug	\$55.00
A00791	Mouse Anti-FES Monoclonal Antibody	40 ug	\$55.00
A00873	Rabbit Anti-KT3 [HRP] Polyclonal Antibody	40 ug	\$55.00
A00820	Mouse Anti-Myoglobin Monoclonal Antibody	40 ug	\$55.00
A00786	Mouse Anti-TYRO3 Monoclonal Antibody	40 ug	\$55.00
A00746	Mouse Anti-CHK1 Monoclonal Antibody	40 ug	\$58.00
A00160	Goat Anti-Mouse IgG (H&L) [HRP] Polyclonal Antibody	1 mg	\$59.50
A00098	Goat Anti-Rabbit IgG (H&L) [HRP] Polyclonal Antibody	1 mg	\$59.50
A00166	Goat Anti-Human IgG (H&L) [HRP] Polyclonal Antibody	1 mg	\$59.50
A00167	Goat Anti-Rat IgG (H&L) [HRP] Polyclonal Antibody	1 mg	\$59.50
A00728	Chicken Anti-Protein A Polyclonal Antibody	40 ug	\$60.00
A00192	Goat Anti-GAPDH [HRP] Polyclonal Antibody	40 ug	\$60.00
A00915	Goat Anti-GAPDH [Biotin] Polyclonal Antibody	40 ug	\$60.00
A00165	Goat Anti-Chicken IgG(H&L) [HRP] Polyclonal Antibody	1 mg	\$63.00
A00193	Rabbit Anti-GSK3β (Phospho-Ser9) Polyclonal Antibody	40 ug	\$69.00
A00702	THETM Anti-β-actin Monoclonal Antibody (Mouse)	40 ug	\$69.00
A00865	THETM Anti-GST Monoclonal Antibody (Mouse)	40 ug	\$69.00
A00337	Rabbit Anti-HDAC5 (Phospho-Ser498) Polyclonal Antibody	40 ug	\$69.00
A00563	Rabbit Anti-CaMKII (Phospho-Thr286) Polyclonal Antibody	40 ug	\$69.00
A00704	THETM Anti-c-Myc-tag Monoclonal Antibody (Mouse)	40 ug	\$69.00
A00443	Rabbit Anti-MKK6 (Phospho-Ser207) Polyclonal Antibody	40 ug	\$69.00
A00385	Rabbit Anti-Akt2 (Phospho-Ser474) Polyclonal Antibody	40 ug	\$69.00
A00402	Rabbit Anti-Tau (Phospho-Ser235) Polyclonal Antibody	40 ug	\$69.00
A00437	Rabbit Anti-ATM (Phospho-Ser1981) Polyclonal Antibody	40 ug	\$69.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00252	Rabbit Anti-STAT4 (Phospho-Tyr693) Polyclonal Antibody	40 ug	\$69.00
A00272	Rabbit Anti-Akt (Phospho-Ser473) Polyclonal Antibody	40 ug	\$69.00
A00250	Rabbit Anti-STAT3 (Phospho-Tyr705) Polyclonal Antibody	40 ug	\$69.00
A00371	Rabbit Anti-IGF-1R (Phospho-Tyr1161) Polyclonal Antibody	40 ug	\$69.00
A00326	Rabbit Anti-Stathmin 1 (Phospho-Ser24) Polyclonal Antibody	40 ug	\$69.00
A00559	Rabbit Anti-eIF2α (Phospho-Ser51) Polyclonal Antibody	40 ug	\$69.00
A00384	Rabbit Anti-c-Abl (Phospho-Tyr412) Polyclonal Antibody	40 ug	\$69.00
A00466	Rabbit Anti-c-Jun (Phospho-Ser63) Polyclonal Antibody	40 ug	\$69.00
A00289	Rabbit Anti-Chk2 (Phospho-Thr68) Polyclonal Antibody	40 ug	\$69.00
A00649	Rabbit Anti-PKR (Phospho-Thr446) Polyclonal Antibody	40 ug	\$69.00
A00399	Rabbit Anti-Tau (Phospho-Ser396) Polyclonal Antibody	40 ug	\$69.00
A00483	Rabbit Anti-p53 (Phospho-Ser15) Polyclonal Antibody	40 ug	\$69.00
A00194	Rabbit Anti-c-Jun (Phospho-Ser73) Polyclonal Antibody	40 ug	\$69.00
A00324	Rabbit Anti-Stathmin 1 (Phospho-Ser37) Polyclonal Antibody	40 ug	\$69.00
A00275	Rabbit Anti-Akt (Phospho-Thr308) Polyclonal Antibody	40 ug	\$69.00
A00338	Rabbit Anti-HDAC4 (Phospho-Ser632) Polyclonal Antibody	40 ug	\$69.00
A00267	Rabbit Anti-CREB (Phospho-Ser133) Polyclonal Antibody	40 ug	\$69.00
A00526	Rabbit Anti-p44/42 MAP Kinase (Phospho-Tyr204) Polyclonal Antibody	40 ug	\$69.00
A00347	Rabbit Anti-STAT1 (Phospho-Ser727) Polyclonal Antibody	40 ug	\$69.00
A00366	Rabbit Anti-VEGFR2 (Phospho-Tyr1175) Polyclonal Antibody	40 ug	\$69.00
A00277	Rabbit Anti-PTEN (Phospho-Ser380, Thr382, 383) Polyclonal Antibody	40 ug	\$69.00
A00296	Rabbit Anti-BAD (Phospho-Ser155) Polyclonal Antibody	40 ug	\$69.00
A00207	Rabbit Anti-PDK1 (Phospho-Ser241) Polyclonal Antibody	40 ug	\$69.00
A00249	Rabbit Anti-STAT1 (Phospho-Tyr701) Polyclonal Antibody	40 ug	\$69.00
A00537	Rabbit Anti-p38 MAPK (Phospho-Tyr182) Polyclonal Antibody	40 ug	\$69.00
A00356	Rabbit Anti-IKKα (Phospho-Thr23) Polyclonal Antibody	40 ug	\$69.00
A00535	Rabbit Anti-p38 MAPK (Phospho-Thr180) Polyclonal Antibody	40 ug	\$69.00
A00376	Rabbit Anti-IGF-1R (Phospho-Tyr1165, 1166) Polyclonal Antibody	40 ug	\$69.00
A00246	Rabbit Anti-Elk-1 (Phospho-Thr417) Polyclonal Antibody	40 ug	\$69.00
A00242	Rabbit Anti-Myc (Phospho-Thr58) Polyclonal Antibody	40 ug	\$69.00
A00505	Rabbit Anti-Stathmin 1 (Phospho-Ser15) Polyclonal Antibody	40 ug	\$69.00
A00504	Rabbit Anti-eIF4E (Phospho-Ser209) Polyclonal Antibody	40 ug	\$69.00
A00369	Rabbit Anti-cofilin (Phospho-Ser3) Polyclonal Antibody	40 ug	\$69.00
A00525	Rabbit Anti-p44/42 MAP Kinase (Phospho-Thr202) Polyclonal Antibody	40 ug	\$69.00
A00659	Rabbit Anti-PLCγ2 (Phospho-Tyr753) Polyclonal Antibody	40 ug	\$69.00
A00236	Rabbit Anti-c-Jun (Phospho-Thr239) Polyclonal Antibody	40 ug	\$69.00
A00339	Rabbit Anti-Histone H3.1 (Phospho-Ser10) Polyclonal Antibody	40 ug	\$69.00
A00379	Rabbit Anti-Paxillin (Phospho-Tyr118) Polyclonal Antibody	40 ug	\$69.00
A00271	Rabbit Anti-EGFR (Phospho-Tyr869) Polyclonal Antibody	40 ug	\$69.00
A00377	Rabbit Anti-Rb (Phospho-Ser780) Polyclonal Antibody	40 ug	\$69.00
A00420	Rabbit Anti-cdc25C (Phospho-Ser216) Polyclonal Antibody	40 ug	\$69.00
A00224	Rabbit Anti-NF-κB p65 (Phospho-Ser536) Polyclonal Antibody	40 ug	\$69.00
A00395	Rabbit Anti-FKHR (Phospho-Ser256) Polyclonal Antibody	40 ug	\$69.00
A00220	Rabbit Anti-NF-κB p65 (Phospho-Ser276) Polyclonal Antibody	40 ug	\$69.00
A00477	Rabbit Anti-MEF2A (Phospho-Thr312) Polyclonal Antibody	40 ug	\$69.00
A00365	Rabbit Anti-Chk2 (Phospho-Ser516) Polyclonal Antibody	40 ug	\$69.00
A00343	Rabbit Anti-HSP27 (Phospho-Ser15) Polyclonal Antibody	40 ug	\$69.00
A00341	Rabbit Anti-ASK1 (Phospho-Ser83) Polyclonal Antibody	40 ug	\$69.00
A00362	Rabbit Anti-JAK2 (Phospho-Tyr221) Polyclonal Antibody	40 ug	\$69.00
A00318	Rabbit Anti-HER2 (Phospho-Tyr1248) Polyclonal Antibody	40 ug	\$69.00
A00541	Rabbit Anti-Connexin 43 (Phospho-Ser367) Polyclonal Antibody	40 ug	\$69.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00405	Rabbit Anti-Tau (Phospho-Thr205) Polyclonal Antibody	40 ug	\$69.00
A00469	Rabbit Anti-NF-κB p100/p52 (Phospho-Ser865) Polyclonal Antibody	40 ug	\$69.00
A00494	Rabbit Anti-DNA PKcs (Phospho-Thr2609) Polyclonal Antibody	40 ug	\$69.00
A00403	Rabbit Anti-Tau (Phospho-Thr181) Polyclonal Antibody	40 ug	\$69.00
A00378	Rabbit Anti-Rb (Phospho-Ser807) Polyclonal Antibody	40 ug	\$69.00
A00291	Rabbit Anti-Ezrin (Phospho-Tyr353) Polyclonal Antibody	40 ug	\$69.00
A00269	Rabbit Anti-ATF-4 (Phospho-Ser245) Polyclonal Antibody	40 ug	\$69.00
A00489	Rabbit Anti-Tau (Phospho-Ser262) Polyclonal Antibody	40 ug	\$69.00
A00290	Rabbit Anti-PTEN (Phospho-Ser370) Polyclonal Antibody	40 ug	\$69.00
A00333	Rabbit Anti-p27Kip1 (Phospho-Thr187) Polyclonal Antibody	40 ug	\$69.00
A00556	Rabbit Anti-CREB (Phospho-Ser129) Polyclonal Antibody	40 ug	\$69.00
A00643	Rabbit Anti-mTOR (Phospho-Ser2448) Polyclonal Antibody	40 ug	\$69.00
A00374	Rabbit Anti-CDC2 (Phospho-Thr161) Polyclonal Antibody	40 ug	\$69.00
A00396	Rabbit Anti-Src (Phospho-Tyr418) Polyclonal Antibody	40 ug	\$69.00
A00665	Rabbit Anti-AMPK1 (Phospho-Thr174) Polyclonal Antibody	40 ug	\$69.00
A00552	Rabbit Anti-Histone H2A.X (Phospho-Ser139) Polyclonal Antibody	40 ug	\$69.00
A00349	Rabbit Anti-IκB-α (Phospho-Tyr42) Polyclonal Antibody	40 ug	\$69.00
A00283	Rabbit Anti-p95/NBS1 (Phospho-Ser343) Polyclonal Antibody	40 ug	\$69.00
A00435	Rabbit Anti-Chk1 (Phospho-Ser345) Polyclonal Antibody	40 ug	\$69.00
A00548	Rabbit Anti-EGFR (Phospho-Tyr1110) Polyclonal Antibody	40 ug	\$69.00
A00570	Rabbit Anti-MARCKS (Phospho-Ser158) Polyclonal Antibody	40 ug	\$69.00
A00367	Rabbit Anti-VEGFR2 (Phospho-Tyr1214) Polyclonal Antibody	40 ug	\$69.00
A00524	Rabbit Anti-CDC2 (Phospho-Tyr15) Polyclonal Antibody	40 ug	\$69.00
A00216	Rabbit Anti-GSK3α (Phospho-Ser21) Polyclonal Antibody	40 ug	\$69.00
A00345	Rabbit Anti-AMPK1/AMPK2 (Phospho-Ser485/491) Polyclonal Antibody	40 ug	\$69.00
A00408	Rabbit Anti-Tau (Phospho-Ser404) Polyclonal Antibody	40 ug	\$69.00
A00342	Rabbit Anti-SEK1/MKK4 (Phospho-Thr261) Polyclonal Antibody	40 ug	\$69.00
A00498	Rabbit Anti-Pyk2 (Phospho-Tyr402) Polyclonal Antibody	40 ug	\$69.00
A00320	Rabbit Anti-EGFR (Phospho-Ser1070) Polyclonal Antibody	40 ug	\$69.00
A00472	Rabbit Anti-JunD (Phospho-Ser255) Polyclonal Antibody	40 ug	\$69.00
A00727	Rabbit Anti-Chk1 (Phospho-Ser296) Polyclonal Antibody	40 ug	\$69.00
A00254	Rabbit Anti-STAT5A (Phospho-Ser780) Polyclonal Antibody	40 ug	\$69.00
A00340	Rabbit Anti-ASK1 (Phospho-Ser966) Polyclonal Antibody	40 ug	\$69.00
A00406	Rabbit Anti-Tau (Phospho-Thr231) Polyclonal Antibody	40 ug	\$69.00
A00496	Rabbit Anti-Bcr (Phospho-Tyr177) Polyclonal Antibody	40 ug	\$69.00
A00232	Rabbit Anti-c-Jun (Phospho-Thr91) Polyclonal Antibody	40 ug	\$69.00
A00654	Rabbit Anti-PKCβ (Phospho-Thr641) Polyclonal Antibody	40 ug	\$69.00
A00361	Rabbit Anti-EGFR (Phospho-Tyr1092) Polyclonal Antibody	40 ug	\$69.00
A00253	Rabbit Anti-STAT5A (Phospho-Tyr694) Polyclonal Antibody	40 ug	\$69.00
A00382	Rabbit Anti-HDAC8 (Phospho-Ser39) Polyclonal Antibody	40 ug	\$69.00
A00295	Rabbit Anti-BAD (Phospho-Ser112) Polyclonal Antibody	40 ug	\$69.00
A00381	Rabbit Anti-Rb (Phospho-Ser795) Polyclonal Antibody	40 ug	\$69.00
A00650	Rabbit Anti-PAK1 (Phospho-Thr212) Polyclonal Antibody	40 ug	\$69.00
A00359	Rabbit Anti-IκB-α (Phospho-Ser32, 36) Polyclonal Antibody	40 ug	\$69.00
A00251	Rabbit Anti-STAT3 (Phospho-Ser727) Polyclonal Antibody	40 ug	\$69.00
A00229	Rabbit Anti-NF-κB p105/p50 (Phospho-Ser893) Polyclonal Antibody	40 ug	\$69.00
A00292	Rabbit Anti-BCL-2 (Phospho-Thr56) Polyclonal Antibody	40 ug	\$69.00
A00648	Rabbit Anti-p62Dok (Phospho-Tyr398) Polyclonal Antibody	40 ug	\$69.00
A00313	Rabbit Anti-HER2 (Phospho-Tyr877) Polyclonal Antibody	40 ug	\$69.00
A00334	Rabbit Anti-Ezrin (Phospho-Thr566) Polyclonal Antibody	40 ug	\$69.00
A00644	Rabbit Anti-Cortactin (Phospho-Tyr466) Polyclonal Antibody	40 ug	\$69.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00579	Rabbit Anti-PKCθ (Phospho-Ser676) Polyclonal Antibody	40 ug	\$69.00
A00668	Rabbit Anti-EGFR (Phospho-Thr693) Polyclonal Antibody	40 ug	\$69.00
A00375	Rabbit Anti-CDK2 (Phospho-Thr160) Polyclonal Antibody	40 ug	\$69.00
A00667	Rabbit Anti-EGFR (Phospho-Thr678) Polyclonal Antibody	40 ug	\$69.00
A00532	Rabbit Anti-SAPK/JNK (Phospho-Thr183) Polyclonal Antibody	40 ug	\$69.00
A00287	Rabbit Anti-FAK (Phospho-Tyr861) Polyclonal Antibody	40 ug	\$69.00
A00439	Rabbit Anti-cdc25A (Phospho-Ser75) Polyclonal Antibody	40 ug	\$69.00
A00330	Rabbit Anti-NF-κB p65 (Phospho-Ser529) Polyclonal Antibody	40 ug	\$69.00
A00484	Rabbit Anti-p53 (Phospho-Ser46) Polyclonal Antibody	40 ug	\$69.00
A00329	Rabbit Anti-4E-BP1 (Phospho-Thr36) Polyclonal Antibody	40 ug	\$69.00
A00285	Rabbit Anti-p73 (Phospho-Tyr99) Polyclonal Antibody	40 ug	\$69.00
A00263	Rabbit Anti-STAT6 (Phospho-Tyr641) Polyclonal Antibody	40 ug	\$69.00
A00394	Rabbit Anti-Caveolin-1(Phospho-Tyr14) Polyclonal Antibody	40 ug	\$69.00
A00240	Rabbit Anti-JunB (Phospho-Ser259) Polyclonal Antibody	40 ug	\$69.00
A00218	Rabbit Anti-PTEN (Phospho-Ser380) Polyclonal Antibody	40 ug	\$69.00
A00549	Rabbit Anti-Merlin (Phospho-Ser518) Polyclonal Antibody	40 ug	\$69.00
A00368	Rabbit Anti-Vav (Phospho-Tyr174) Polyclonal Antibody	40 ug	\$69.00
A00502	Rabbit Anti-EGFR (Phospho-Tyr1197) Polyclonal Antibody	40 ug	\$69.00
A00501	Rabbit Anti-EGFR (Phospho-Tyr1172) Polyclonal Antibody	40 ug	\$69.00
A00237	Rabbit Anti-c-Jun (Phospho-Ser243) Polyclonal Antibody	40 ug	\$69.00
A00387	Rabbit Anti-FAK (Phospho-Tyr925) Polyclonal Antibody	40 ug	\$69.00
A00546	Rabbit Anti-HSF1 (Phospho-Ser303) Polyclonal Antibody	40 ug	\$69.00
A00568	Rabbit Anti-Gab1 (Phospho-Tyr627) Polyclonal Antibody	40 ug	\$69.00
A00475	Rabbit Anti-ATF-2 (Phospho-Thr73/55) Polyclonal Antibody	40 ug	\$69.00
A00321	Rabbit Anti-Met (Phospho-Tyr1234) Polyclonal Antibody	40 ug	\$69.00
A00566	Rabbit Anti-BIM (Phospho-Ser65/69) Polyclonal Antibody	40 ug	\$69.00
A00543	Rabbit Anti-NF-κB p65 (Phospho-Ser311) Polyclonal Antibody	40 ug	\$69.00
A00497	Rabbit Anti-MEK1 (Phospho-Ser217) Polyclonal Antibody	40 ug	\$69.00
A01024	Rabbit Anti-Lck (Phospho-Tyr393) Polyclonal Antibody	40 ug	\$69.00
A00297	Rabbit Anti-Estrogen Receptor-α (Phospho-Ser106) Polyclonal Antibody	40 ug	\$69.00
A00657	Rabbit Anti-PKCθ (Phospho-Ser695) polyclonal antibody	40 ug	\$69.00
A00471	Rabbit Anti-c-Jun (Phospho-Tyr170) Polyclonal Antibody	40 ug	\$69.00
A01026	Rabbit Anti-APP (Phospho-Thr668) Polyclonal Antibody	40 ug	\$69.00
A00470	Rabbit Anti-NF-κB p105/p50 (Phospho-Ser337) Polyclonal Antibody	40 ug	\$69.00
A00495	Rabbit Anti-SEK1/MKK4 (Phospho-Ser80) Polyclonal Antibody	40 ug	\$69.00
A01027	Rabbit Anti-IRS-1 (Phospho-Ser639) Polyclonal Antibody	40 ug	\$69.00
A00360	Rabbit Anti-JAK2 (Phospho-Tyr1007) Polyclonal Antibody	40 ug	\$69.00
A00540	Rabbit Anti-RelB (Phospho-Ser552) Polyclonal Antibody	40 ug	\$69.00
A00230	Rabbit Anti-NF-κB p105/p50 (Phospho-Ser907) Polyclonal Antibody	40 ug	\$69.00
A01028	Rabbit Anti-GluR1 (Phospho-Ser849) Polyclonal Antibody	40 ug	\$69.00
A00561	Rabbit Anti-p70 S6 Kinase (Phospho-Ser424) Polyclonal Antibody	40 ug	\$69.00
A00294	Rabbit Anti-BCL-XL (Phospho-Ser62) Polyclonal Antibody	40 ug	\$69.00
A00315	Rabbit Anti-S6 Ribosomal Protein (Phospho-Ser235) Polyclonal Antibody	40 ug	\$69.00
A00468	Rabbit Anti-NF-κB p65 (Phospho-Thr254) Polyclonal Antibody	40 ug	\$69.00
A00539	Rabbit Anti-p70 S6 Kinase (Phospho-Thr421) Polyclonal Antibody	40 ug	\$69.00
A00492	Rabbit Anti-Zap-70 (Phospho-Tyr493) Polyclonal Antibody	40 ug	\$69.00
A01029	Rabbit Anti-p56Dok-2 (Phospho-Tyr299) Polyclonal Antibody	40 ug	\$69.00
A00293	Rabbit Anti-BCL-2 (Phospho-Ser70) Polyclonal Antibody	40 ug	\$69.00
A00358	Rabbit Anti-Src (Phospho-Tyr529) Polyclonal Antibody	40 ug	\$69.00
A00314	Rabbit Anti-HER2 (Phospho-Tyr1221, 1222) Polyclonal Antibody	40 ug	\$69.00
A00336	Rabbit Anti-SMC1 (Phospho-Ser957) Polyclonal Antibody	40 ug	\$69.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00491	Rabbit Anti-Zap-70 (Phospho-Tyr319) Polyclonal Antibody	40 ug	\$69.00
A01030	Rabbit Anti-PKR (Phospho-Thr451) Polyclonal Antibody	40 ug	\$69.00
A00335	Rabbit Anti-Paxillin (Phospho-Tyr31) Polyclonal Antibody	40 ug	\$69.00
A00490	Rabbit Anti-BRCA1 (Phospho-Ser1524) Polyclonal Antibody	40 ug	\$69.00
A01031	Rabbit Anti-GSK3β (Phospho-Tyr216) Polyclonal Antibody	40 ug	\$69.00
A00647	Rabbit Anti-p62Dok (phospho-Tyr362) Polyclonal Antibody	40 ug	\$69.00
A00924	Rabbit Anti-CARM1 (Phospho-Ser228) Polyclonal Antibody	40 ug	\$69.00
A00248	Rabbit Anti-GATA1 (Phospho-Ser310) Polyclonal Antibody	40 ug	\$69.00
A00312	Rabbit Anti-Progesterone Receptor (Phospho-Ser190) Polyclonal Antibody	40 ug	\$69.00
A00697	Rabbit Anti-APP (Phospho-Tyr757) Polyclonal Antibody	40 ug	\$69.00
A00226	Rabbit Anti-NF-κB p100/p52 (Phospho-Ser869) Polyclonal Antibody	40 ug	\$69.00
A01032	Rabbit Anti-IκB-β (Phospho-Ser23) Polyclonal Antibody	40 ug	\$69.00
A00923	Rabbit Anti-α-Catenin (Phospho-Ser641) Polyclonal Antibody	40 ug	\$69.00
A00247	Rabbit Anti-GATA1 (Phospho-Ser142) Polyclonal Antibody	40 ug	\$69.00
A00969	Rabbit Anti-SHP-1 (Phospho-Tyr536) Polyclonal Antibody	40 ug	\$69.00
A00534	Rabbit Anti-NF-κB p105/p50 (Phospho-Ser932) Polyclonal Antibody	40 ug	\$69.00
A00487	Rabbit Anti-p90RSK (Phospho-Thr348) Polyclonal Antibody	40 ug	\$69.00
A01033	Rabbit Anti-Keratin 18 (Phospho-Ser33) Polyclonal Antibody	40 ug	\$69.00
A00354	Rabbit Anti-FKHRL1 (Phospho-Ser253) Polyclonal Antibody	40 ug	\$69.00
A00922	Rabbit Anti-NFκB-p105/p50 (Phospho-Ser927) Polyclonal Antibody	40 ug	\$69.00
A00332	Rabbit Anti-VASP (Phospho-Ser157) Polyclonal Antibody	40 ug	\$69.00
A00486	Rabbit Anti-NMDAR1 (Phospho-Ser896) Polyclonal Antibody	40 ug	\$69.00
A01034	Rabbit Anti-MAPKAPK-2 (Phospho-Thr334) Polyclonal Antibody	40 ug	\$69.00
A00288	Rabbit Anti-Integrin β3 (Phospho-Tyr773) Polyclonal Antibody	40 ug	\$69.00
A00353	Rabbit Anti-VASP (Phospho-Ser238) Polyclonal Antibody	40 ug	\$69.00
A00440	Rabbit Anti-MKK3 (Phospho-Ser189) Polyclonal Antibody	40 ug	\$69.00
A00921	Rabbit Anti-HSP90B (Phospho-Ser254) Polyclonal Antibody	40 ug	\$69.00
A00245	Rabbit Anti-Elk-1 (Phospho-Ser389) Polyclonal Antibody	40 ug	\$69.00
A00485	Rabbit Anti-p53 (Phospho-Ser315) Polyclonal Antibody	40 ug	\$69.00
A00223	Rabbit Anti-NF-κB p65 (Phospho-Ser468) Polyclonal Antibody	40 ug	\$69.00
A01035	Rabbit Anti-Shc1 (Phospho-Tyr427) Polyclonal Antibody	40 ug	\$69.00
A00554	Rabbit Anti-p70 S6 Kinase (Phospho-Ser411) Polyclonal Antibody	40 ug	\$69.00
A00920	Rabbit Anti-ADD1 (Phospho-Ser726) Polyclonal Antibody	40 ug	\$69.00
A00244	Rabbit Anti-Myc (Phospho-Ser373) Polyclonal Antibody	40 ug	\$69.00
A00373	Rabbit Anti-FKHR (Phospho-Ser319) Polyclonal Antibody	40 ug	\$69.00
A00965	Rabbit Anti-Akt (Phospho-Ser473) Polyclonal Antibody	40 ug	\$69.00
A01036	Rabbit Anti-Smad1 (Phospho-Ser465) Polyclonal Antibody	40 ug	\$69.00
A00351	Rabbit Anti-MEK1 (Phospho-Ser221) Polyclonal Antibody	40 ug	\$69.00
A00894	Rabbit Anti-Tau (Phospho-Ser199) Polyclonal Antibody	40 ug	\$69.00
A00919	Rabbit Anti-14-3-3 ζ (Phospho-Ser58) Polyclonal Antibody	40 ug	\$69.00
A00507	Rabbit Anti-Met (Phospho-Tyr1349) Polyclonal Antibody	40 ug	\$69.00
A00243	Rabbit Anti-Myc (Phospho-Thr358) Polyclonal Antibody	40 ug	\$69.00
A00576	Rabbit Anti-PKCδ (Phospho-Ser645) Polyclonal Antibody	40 ug	\$69.00
A00372	Rabbit Anti-AFX (Phospho-Ser197) Polyclonal Antibody	40 ug	\$69.00
A00195	Rabbit Anti-Elk-1 (Phospho-Ser383) Polyclonal Antibody	40 ug	\$69.00
A00530	Rabbit Anti-HSP27 (Phospho-Ser82) Polyclonal Antibody	40 ug	\$69.00
A00221	Rabbit Anti-NF-κB p65 (Phospho-Thr435) Polyclonal Antibody	40 ug	\$69.00
A01037	Rabbit Anti-Smad2 (Phospho-Ser467) Polyclonal Antibody	40 ug	\$69.00
A00506	Rabbit Anti-IRS-1 (Phospho-Ser307) Polyclonal Antibody	40 ug	\$69.00
A00328	Rabbit Anti-4E-BP1 (Phospho-Thr45) Polyclonal Antibody	40 ug	\$69.00
A00481	Rabbit Anti-PKD/PKCμ (Phospho-Ser738) Polyclonal Antibody	40 ug	\$69.00

MiniFresh® Antibodies:

Small packages to help you to reduce your research costs.



Cat. No.	Name	Size	Price
A00551	Rabbit Anti-Synapsin (Phospho-Ser9) Polyclonal Antibody	40 ug	\$69.00
A00241	Rabbit Anti-ATF-2 (Phospho-Thr71/53) Polyclonal Antibody	40 ug	\$69.00
A00663	Rabbit Anti-Integrin β3 (phospho-Tyr785) Polyclonal Antibody	40 ug	\$69.00
A00370	Rabbit Anti-VEGFR2 (Phospho-Tyr951) Polyclonal Antibody	40 ug	\$69.00
A00528	Rabbit Anti-HSP27 (Phospho-Ser78) Polyclonal Antibody	40 ug	\$69.00
A00480	Rabbit Anti-MyoD (Phospho-Ser200) Polyclonal Antibody	40 ug	\$69.00
A00479	Rabbit Anti-Estrogen Receptor-α (Phospho-Ser104) Polyclonal Antibody	40 ug	\$69.00
A00503	Rabbit Anti-IRS-1 (Phospho-Ser636) Polyclonal Antibody	40 ug	\$69.00
A00239	Rabbit Anti-JunB (Phospho-Ser79) Polyclonal Antibody	40 ug	\$69.00
A00571	Rabbit Anti-MEK1 (Phospho-Thr291) Polyclonal Antibody	40 ug	\$69.00
A00898	Rabbit Anti-Tau (Phospho-Ser400) Polyclonal Antibody	40 ug	\$69.00
A00478	Rabbit Anti-STAT6 (Phospho-Thr645) Polyclonal Antibody	40 ug	\$69.00
A00217	Rabbit Anti-MEK2 (Phospho-Thr394) Polyclonal Antibody	40 ug	\$69.00
A00523	Rabbit Anti-BRCA1 (Phospho-Ser1423) Polyclonal Antibody	40 ug	\$69.00
A00417	Rabbit Anti-β-Catenin (Phospho-Thr41, Ser45) Polyclonal Antibody	40 ug	\$69.00
A00476	Rabbit Anti-ATF-2 (Phospho-Ser112/94) Polyclonal Antibody	40 ug	\$69.00
A00215	Rabbit Anti-Raf-1 (Phospho-Ser259) Polyclonal Antibody	40 ug	\$69.00
A00913	Rabbit Anti-PKCα/β2 (Phospho-Thr638/641) Polyclonal Antibody	40 ug	\$69.00
A00500	Rabbit Anti-β-Catenin (Phospho-Ser37) Polyclonal Antibody	40 ug	\$69.00
A01020	Rabbit Anti-ICAM-1 (Phospho-Tyr512) Polyclonal Antibody	40 ug	\$69.00
A00499	Rabbit Anti-β-Catenin (Phospho-Ser33) Polyclonal Antibody	40 ug	\$69.00
A01021	Rabbit Anti-PKD/PKCμ (Phospho-Ser910) Polyclonal Antibody	40 ug	\$69.00
A00299	Rabbit Anti-Estrogen Receptor-α (Phospho-Ser167) Polyclonal Antibody	40 ug	\$69.00
A00364	Rabbit Anti-G3BP-1 (Phospho-Ser232) Polyclonal Antibody	40 ug	\$69.00
A00474	Rabbit Anti-ATF-2 (Phospho-Thr69/51) Polyclonal Antibody	40 ug	\$69.00
A00896	Rabbit Anti-Tau (Phospho-Thr217) Polyclonal Antibody	40 ug	\$69.00
A00234	Rabbit Anti-c-Jun (Phospho-Thr93) Polyclonal Antibody	40 ug	\$69.00
A01022	Rabbit Anti-PLCγ1 (Phospho-Tyr783) Polyclonal Antibody	40 ug	\$69.00
A00298	Rabbit Anti-Estrogen Receptor-α (Phospho-Ser118) Polyclonal Antibody	40 ug	\$69.00
A00363	Rabbit Anti-JAK1 (Phospho-Tyr1022) Polyclonal Antibody	40 ug	\$69.00
A00730	THETM Anti-β-Actin [HRP] Monoclonal Antibody (Mouse)	40 ug	\$74.00
A00866	THETM Anti-GST [HRP] Monoclonal Antibody (Mouse)	40 ug	\$74.00
A00864	THETM Anti-c-Myc [Biotin] Monoclonal Antibody (Mouse)	40 ug	\$74.00
A00867	THETM Anti-GST-tag [Biotin] Monoclonal Antibody (Mouse)	40 ug	\$74.00
A00863	THETM Anti-c-Myc [HRP] Monoclonal Antibody (Mouse)	40 ug	\$74.00