

CERTIFICATE OF ANALYSIS

Product information

Product name	Goat anti Rabbit IgG Magbeads
Cat. No.	L00328
Lot No.	R314904
Storage condition	Store at 4°C

QC results

Test Items	Specifications	Results
Binding capacity	$\geq 150\mu\text{g/ml}$	Pass
Stability	Binding capacity keep above 20mg/ml after incubate at 37°C for 1 week.	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: *Richard* Date: 04/17/2009
Quality Assurance Department