

CERTIFICATE OF ANALYSIS

Product Information

Product name:	Taq DNA polymerase
Cat. No.:	E00007
Lot No.:	R204907
Storage condition:	Store at -20°C

QC Results

Test Items	Specifications	Results
Activity	5U/ul	Pass
Amplification capability	Up to 8000bp DNA fragment.	Pass
3.8Kb fragment amplification	Amplifying 3.8Kb beta-globin gene from human genomic DNA	Pass
Exo - Nuclease contamination	No contamination	Pass
<i>E.coli</i> Genomic DNA contamination	No contamination	Pass
Plasmid DNA contamination	No contamination	Pass
Mycoplasma contamination	No contamination	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:



Date: 07/28/2009

Quality Assurance Department