

## CERTIFICATE OF ANALYSIS

### Product Information

---

Product name:	Taq DNA polymerase(2500U)
Cat. No.:	E00007
Lot No.:	R219906
Storage condition:	Store at -20°C

---

### QC Results

---

Test Items	Specifications	Results
Activity	5U/ul	Pass
Amplification capability	Up to 8000bp DNA fragment.	Pass
Endo - Nuclease contamination	No contamination	Pass
Exo - Nuclease contamination	No contamination	Pass
<i>E.coli</i> Genomic DNA contamination	No contamination	Pass
Plasmid DNA contamination	No contamination	Pass
Mycoplasma contamination	No contamination	Pass

---

### Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  
Quality Assurance Department



Date: 06/30/2009