


Datasheet Update: 03262009

## pRNA-H1.1/Hygro

Cat. No. SD1204

**Description:** pRNA-H1.1/Hygro is a GenScript siRNA expression vector. It is designed for mammalian transfection. It carries a Hygromycin resistance gene which can be used for establishing stable cell line. It uses H1 promoter for siRNA expression.



**Polylinker:** 4819-4842  
**H1 Promoter:** 4719-4818  
**SV40 Promoter:** 597-942  
**Hygromycin:** 965-1990  
**pUC ori:** 2660-3300  
**Ampicillin:** 3448-4308

**Forward Sequencing Primer:**

[DA0013: pRNA-H1 Forward](#)  
(TAATACGACTCACTATAGGG)

**Reverse Sequencing Primer:**

[DA0012: pRNA Reverse](#)  
(TAGAAGGCACAGTCGAGG)

T7 Primer → GGATCCGAGCTCGGTACCAAGCTTAAGTTAAACCGCTGATCAGCCTCGACTGTGCCTTCTAGT  
 BamHI Asp718 I Hind III Kpn I BGH reverse primer site