




siRNA expression vector pRNA-H1.1/Shuttle (Cat. No.: SD1206)

Description: pRNA-H1.1/Shuttle is a GenScript adenoviral siRNA shuttle vector. It is compatible with Clontech Adeno-X™ Expression System 1. It uses H1 promoter for siRNA expression.



Polylinker: 4753-4776
H1 Promoter: 4653-4752
SV40 Promoter: 674-1019
Neomycin: 1060-1854
pUC ori: 2568-3208
Ampicillin: 3356-4216

Forward Sequencing Primer:
[DA0013: pRNA-H1 Forward](#)
 (TAATACGACTCACTATAGGG)
Reverse Sequencing Primer:
[DA0012: pRNA Reverse](#)
 (TAGAAGGCACAGTCGAGG)

