


Datasheet Update: 03272009

pRNA-H1.1/Zeo

Cat. No. SD1208

Description: pRNA-H1.1/Zeo is a GenScript siRNA expression vector. It is designed for mammalian transfection. It carries a Zeomycin resistance gene which can be used for establishing stable cell line. It uses H1 promoter for siRNA expression.



Polylinker: 375-398
H1 Promoter: 275-374
SV40 Promoter: 1001-1347
Zeomycin: 1435-1809
pUC ori: 2479-3122
Ampicillin: 3270-4130

Forward Sequencing Primer:
[DA0013: pRNA-H1 Forward](#)
 (TAATACGACTCACTATAGGG)

Reverse Sequencing Primer:
[DA0012: pRNA Reverse](#)
 (TAGAAGGCACAGTCGAGG)