


## pRNA-H1.1/Adeno

Cat. No. SD1209

**Description:** GenScript pRNA-H1.1/Adeno is an adenoviral siRNA shuttle vector. It is compatible with Stratagene AdEasy Adenoviral Vector System. The vector is designed for transfection in mammalian cells. siRNA insert can be easily cloned into polylinker region under H1 promoter control. The vector also contains a Kanamycin resistance gene for transformant selection. LITR and RITR regions are available in this vector as a function for recombinant viral DNA replication.



**Polylinker:** 506-529

**H1 Promoter:** 406-505

**pBR322 origin:** 4028-4695

**Kanamycin:** 5504-6295

**Forward Sequencing Primer:**

[DA0013\\_pRNA-H1 Forward](#)

(TAATACGACTCACTATAGGG)

**Reverse Sequencing Primer:**

[DA0014\\_pShuttle-H1.1 Reverse](#)

(CAAACTACATAAGACCCCCAC)