

CERTIFICATE OF ANALYSIS

Product Information

Product name:	2xTaq Master Mix
Cat. No.:	E00019
Lot No.:	R201909
Storage condition:	Store at -20°C

QC Result

Test Items	Specifications	Results
Activity	Prepared 2xTaq Master Mix has the same effect as the Taq DNA polymerase got from the USA.	Pass
Amplification capability	Up to 8000bp DNA fragment.	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 09/04/2009
Quality Assurance Department