

CERTIFICATE OF ANALYSIS

Product Information

Product name	Protein L Resin
Cat. No.	L00239
Lot No.	R301911
Storage condition	Store at 4°C

QC Result

Test Item	Specification	Result
Binding capacity	$\geq 15\text{mg/ml}$	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: 
Quality Assurance Department

Date: 11/06/2009