

CERTIFICATE OF ANALYSIS

Product Information

Product name:	2xTaq Master Mix
Cat. No.:	E00019
Lot No.:	R20011001
Storage condition:	Store at -20°C

QC Results

Test Items	Specifications	Results
Activity	Prepared 2xTaq Master Mix has the same effect as the Taq DNA polymerase got from the USA.	Pass
Amplification capability	Up to 8000bp DNA fragment.	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: *Catherine* Date: 01/15/2010
Quality Assurance Department