

CERTIFICATE OF ANALYSIS

Product Information

Product name:	Green Taq DNA polymerase
Cat. No.:	L00342
Lot No.:	R20051001
Storage condition:	Store at -20°C

QC Results

Test Items	Specifications	Results
Activity	5U/ul	Pass
Amplification capability	Up to 8000bp DNA fragment.	Pass
Stability	No activity lost after incubate at 37°C for 1 week.	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:



Quality Assurance Department

Date: 01/15/2010