


## pRNA-H1.1/Neo/siFLuc

Cat. No. SD1503

**Description:** pRNA-H1.1/Neo/siFLuc is a GenScript siRNA expression vector. It is designed for mammalian transfection. It uses H1 promoter for siRNA expression. This vector contains siRNA construct for Firefly luciferase, and can be used as a positive control. It can also be used as a siRNA vector using BamH I and Hind III for siRNA insertion.



**Polylinker:** 4700-4769  
**H1 Promoter:** 4600-4699  
**SV40 Promoter:** 647-992  
**Neomycin:** 1033-1827  
**pUC ori:** 2541-3181  
**Ampicillin:** 3329-4189

**Forward Sequencing Primer:**  
 DA0013: pRNA-H1 Forward  
 (TAATACGACTCACTATAGGG)

**Reverse Sequencing Primer:**  
 DA0012: pRNA Reverse  
 (TAGAAGGCACAGTCGAGG)