


siRNA expression vector pRNA-H1.1/Hygro/siFLuc (Catalogue No.: SD1504)

Description: pRNA-H1.1/Hygro/siFLuc is a GenScript siRNA expression vector. It is designed for mammalian transfection. It uses H1 promoter for siRNA expression. This vector contains siRNA construct for Firefly luciferase, and can be used as a positive control. It can also be used as a siRNA vector using BamH I and Hind III for siRNA insertion.



Polylinker: 4819-4888
H1 Promoter: 4719-4818
SV40 Promoter: 597-942
Hygromycin: 965-1990
pUC ori: 2660-3300
Ampicillin: 3448-4308

Forward Sequencing Primer:
[DA0013: pRNA-H1 Forward](#)
 (TAATACGACTCACTATAGGG)
Reverse Sequencing Primer:
[DA0012: pRNA Reverse](#)
 (TAGAAGGCACAGTCGAGG)

