

CERTIFICATE OF ANALYSIS

Product Information

Product name:	Taq DNA polymerase
Cat. No.:	E00007
Lot No.:	R20011002
Storage condition:	Store at -20°C

QC Results

Test Items	Specifications	Results
Activity	5U/ μ L	PASS
Amplification capability	Up to 8000bp DNA fragment.	PASS
Exo - Nuclease contamination	No contamination	PASS
<i>E.coli</i> Genomic DNA contamination	No contamination	PASS
Plasmid DNA contamination	No contamination	PASS
Mycoplasma contamination	No contamination	PASS

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: 
Quality Assurance Department

Date: 02/09/2010