


Datasheet Update: 10102010

pRNA-H1.3/Neo

Cat. No. SD1552

Description: pRNA-H1.3/Neo is a GenScript siRNA expression vector. It is designed for mammalian transfection. It uses an enhanced H1 promoter (enhanced by a CMV enhancer just before H1 promoter) for siRNA expression. The siRNA insert can be cloned between BamH I and Hind III sites.



Polylinker:	5142-5211
CMV IE:	4559-5042
siFluc:	5148-5205
H1 Promoter:	5049-5143
SV40 Promoter:	647-992
Neomycin:	1033-1827
pUC ori:	2541-3181
Ampicillin:	3329-4189

Forward Sequencing Primer:
[DA0024: pRNA-H1.3 Forward](#)
(GACGTCAATGGGAGTTTGT)

Reverse Sequencing Primer:
[DA0012: pRNA Reverse](#)
(TAGAAGGCACAGTCGAGG)