

## CERTIFICATE OF ANALYSIS

### Product information

Product name	NI-IDA Resin
Cat. No.	L00223
Lot No.	R30011002
Storage condition	Store at 4°C

### QC result

Test Item	Specification	Result
Binding capacity	$\geq 20\text{mg/ml}$	Pass

### Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:   
Quality Assurance Department

Date: 02/08/2010