

CERTIFICATE OF ANALYSIS

Product information

Product name	High Affinity GST Resin
Cat. No.	L00206
Lot No.	R30121002
Storage condition	Store at +4°C to +8°C

QC results

Test Items	Specifications	Results
Binding capacity	≥ 20 mg/ml	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: 
Quality Assurance Department

Date: 02/25/2010