


Positive MS Spectrum



Sample Information

User : Grace  
Date Acquired : 12/26/2009 15:51:58  
Sample Name : Peptide YY (PYY), human  
Cat No. : RP10355  
Lot No. : 85691001121409LQ  
Injection Volum : 1.000µL  
Theoretical MW : 4309.75  
Observed MW : 4311.80

Interface : ESI  
Nebulizing Gas Flow : 4.21 L/min  
CDL Temp : 250°C  
CDL Volt : -25.0V  
Block Temp : 200°C

Prerod Bias : +4.5kv  
Detector : 1.5kv  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH