


**siRNA expression vector pRNA-H1.1/Hygro/siRLuc (Catalogue No.: SD1604)**

**Description:** pRNA-H1.1/Hygro/siRLuc is a GenScript siRNA expression vector. It is designed for mammalian transfection. It uses H1 promoter for siRNA expression. This vector contains siRNA construct for Renilla luciferase, and can be used as a positive control. It can also be used as a siRNA vector using BamH I and Hind III for siRNA insertion.



**Polylinker:** 4819-4888  
**H1 Promoter:** 4719-4818  
**SV40 Promoter:** 597-942  
**Hygromycin:** 965-1990  
**pUC ori:** 2660-3300  
**Ampicillin:** 3448-4308

**Forward Sequencing Primer:**  
[DA0013: pRNA-H1 Forward](#)  
 (TAATACGACTCACTATAGGG)  
**Reverse Sequencing Primer:**  
[DA0012: pRNA Reverse](#)  
 (TAGAAGGCACAGTCGAGG)

