

## CERTIFICATE OF ANALYSIS

### Product information

Product name	Anti-His mAb MagBeads
Cat. No.	L00275
Lot No.	R30021003
Storage condition	Store at +4 to +8°C

### QC result

Test Item	Specification	Result
Binding capacity	$\geq 300$ ug/ml	PASS

### Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 03/02/2010  
Quality Assurance Department