

CERTIFICATE OF ANALYSIS

Product information

Product name	Protein L Resin
Cat. No.	L00239
Lot No.	R30011003
Storage condition	Store at +4 to +8°C

QC results

Test Items	Specifications	Results
Binding capacity	≥ 15 mg/ml	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: 
Quality Assurance Department

Date: 03/02/2010