


## pRNA-H1.1/Neo/sicGFP

Cat. No. SD1702

**Description:** pRNA-H1.1/Neo/sicGFP is a GenScript siRNA expression vector. It is designed for mammalian transfection. It uses H1 promoter for siRNA expression. This vector contains siRNA construct for cGFP (coral green fluorescence protein), and can be used as a positive control. It can also be used as a siRNA vector using BamH I and Hind III for siRNA insertion.



**Polylinker:** 4703-4777

**H1 Promoter:** 4600-4699

**SV40 Promoter:** 647-992

**Neomycin:** 1033-1827

**pUC ori:** 2541-3181

**Ampicillin:** 3329-4189

**Forward Sequencing Primer:**

DA0013: pRNA-H1 Forward  
(TAATACGACTCACTATAGGG)

**Reverse Sequencing Primer:**

DA0012: pRNA Reverse  
(TAGAAGGCACAGTCGAGG)