

CERTIFICATE OF ANALYSIS

Product information

Product name	High-Affinity Antibody Purification Kit
Cat. No.	L00404
Lot No.	R30341003
Storage condition	Store at 4°C

QC results

Test Items	Specifications	Results
Binding capacity	$\geq 3\text{mg/ml}$	Pass

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 03/26/2010
Quality Assurance Department