

CERTIFICATE OF ANALYSIS

Product information

Product name	Mouse Anti-His mAb MagBeads
Cat. No.	L00275
Lot No.	R30121004
Storage condition	Store at +2 to +8°C

QC result

Test Item	Specification	Result
Binding capacity	$\geq 300\text{ug/ml}$	PASS

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 04/08/2010
Quality Assurance Department