

## CERTIFICATE OF ANALYSIS

### Product information

Product name	Protein L Resin
Order No.	L00239
Lot. No.	R30241004
Storage condition	Store at +4 to +8°C

### QC results

Test Items	Specifications	Results
Binding capacity	$\geq 15\text{mg/mL}$	Pass

### Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 04/14/2010  
Quality Assurance Department