

CERTIFICATE OF ANALYSIS

Product Information

Product name	Tau Antibody (Phospho-Thr ²³¹), pAb, Rabbit
Cat. No.	A00406
Lot No.	R50091005
Storage condition	Store at -20°C

QC Result

Test Item	Specification	Result
Concentration(IgG)	1±0.1mg/ml by LOWRY	1.02mg/ml
Purity(IgG)	>90% by SDS-PAGE	Qualified
ELISA	Dilution range: ≥1: 4000	1: 64000
Western Bolt	Dilution range: ≥1: 500 Predicted MW of protein: 48KD	1: 1000 Detected MW of protein: 48KD
Immunohistochemistry	Dilution range: ≥1: 50 Predicted Location of protein: cytoplasmic	1: 100 Detected Location of protein: cytoplasmic Analysis of paraffin-embedded rat hippocampal region tissue

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 05/04/2010
Quality Assurance Department