

CERTIFICATE OF ANALYSIS

Product information

Product Name	MBP-tag Antibody, mAb, Mouse
Product Cat. No.	A00190
Lot No.	10E000140
Quantity	40 µg
Presentation	Lyophilized powder from antibody solution with PBS, pH 7.4, containing 0.02% sodium azide
Conjugation	Unconjugated
Storage	The product is stable in lyophilized form for at least two years if stored at -20°C. The reconstituted antibody can be stored for up to three months at 2-8°C or for up to two years at -20°C or below. Avoid repeated freeze and thaw cycles.
Expiration date	05/25/2013

Quality Control

Test Items	Specifications	Results
Reactivity	ELISA titer $\geq 1:20,000$	Passed
	Specific to MBP fusion protein by Western blot	Passed
Purification	Protein A affinity column	Qualified
Purity	$\geq 85\%$ by SDS-PAGE	Qualified

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: 
Production Supervisor Date: 05/26/2010

Certified by: 
Quality Assurance Manager Date: 05/26/2010