

## PRODUCT INFORMATION

Version: 07/16/2010

### dUTP (100 mM) $\geq$ 99% HPLC

**Cat. No.:** C01717

**Name:** dUTP (100 mM)  $\geq$ 99% HPLC

**Size:** 1 ml, 10 ml, 50 ml and 250 ml

**Cas:** 102814-08-4

**Synonyms:** 2'-Deoxyuridine-5'-Triphosphate, Sodium Salt

**Molecular Formula:**  $C_9H_{12}N_2Na_3O_{14}P_3$

**Molecular Weight:** 534

**Storage:** This product can remain stable for 3 years if stored at  $-20^{\circ}C$ . Always avoid freeze-thaw cycles or exposure to frequent temperature changes.