

CERTIFICATE OF ANALYSIS

Product information

Product Name	Biotin-Protein L
Product Cat. No.	M00097
Lot No.	10I000186
Quantity	1 mg
Presentation	Lyophilized powder from 2ml 0.5 mg/ml Biotin-Protein L solution with PBS, pH 7.4, containing 10 mg/ml BSA and 0.02% sodium azide.
Conjugation	Biotin conjugated
Storage	The product is stable in lyophilized form for at least one year if stored at -20°C. The reconstituted product can be stored for up to three months at 2-8°C or for up to one year at -20°C or below. Avoid repeated freeze and thaw cycles.
Expiration date	9/16/2011

Quality Control

Test Items	Specifications	Results
Reactivity	ELISA titer $\geq 1:20,000$	Passed
	Specific to human IgG by Western blot	Passed

Caution:

For laboratory or further manufacturing use only. Not intended for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: 
Quality Assurance Manager

Date: 9/17/2010