


MS Spectrum

Positive



Sample Information Item
 Date and Time : 2005-3-7 14:46:11
 User : xie
 Sample : RP10094 NV-41-NH2
 Inj. Volume : 1.000
 MW : 4869.55

Cat. No.: RP10094
 Lot No.: PJ00813505

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp: 250°C
 Block Temp: 200°C

Probe bias: +4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. conc: 50%H2O/50%ACN