


Positive MS Spectrum



Sample Information

User : Grace
Date Acquired : 11/26/2009 12:02:30
Sample Name : Pep-1(uncapped)
Cat No. : RP10130
Lot No. : 8369300111909BF
Injection Volum : 1.000µL
Theoretical MW : 2848.30
Observed MW : 2848.05

Interface : ESI
Nebulizing Gas Flow : 1.50 L/min
CDL Temp : 250°C
CDL Volt : -25.0V
Block Temp : 200°C

Prerod Bias : +4.5kv
Detector : 1.5kv
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH