


Positive

### MS Spectrum

Intensity



**Sample Information**

User : Grace  
Date Acquired : 05/23/2009 15:59:19  
Sample Name : MBP  
Cat No. : RP10314  
Lot No. : CC03905  
Injection Volum : 1.000µL  
Theoretical MW : 1390.50  
Observed MW : 1390.00

Interface : ESI  
Nebulizing Gas Flow : 1.50 L/min  
CDL Temp : 250°C  
CDL Volt : -25.0V  
Block Temp : 200°C

Prerod Bias : +4.5kv  
Detector : 1.5kv  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH