


Positive MS Spectrum



Sample Information

User : Grace  
Date Acquired : 11/25/2009 11:46:33  
Sample Name : PACAP (6-38), human, ovine, rat  
Cat No. : RP10338  
Lot No. : 83177002111309ZW  
Injection Volum : 1.000µL  
Theoretical MW : 4024.76  
Observed MW : 4024.00

Interface : ESI  
Nebulizing Gas Flow : 1.50 L/min  
CDL Temp : 250°C  
CDL Volt : -25.0V  
Block Temp : 200°C

Prerod Bias : +4.5kv  
Detector : 1.5kv  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH