


MS Spectrum

Positive



Sample Information Item

Date and Time : 2005-2-2 17:54:55
 User : FJS
 Sample : RP10339(HK-38-NH2)
 Inj. Volume : 1.000
 MW : 4548.38

Cat. No.: RP10339

Lot No.: PJ00613505

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp: 250°C
 Block Temp: 200°C

Probe bias: +4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. conc: 50%H2O/50%ACN