


MS Spectrum

Positive
100



Sample Information Item
Date and Time : 2005-4-5 11:53:59
User : xie
Sample : RP10408 CC-6(0)
Inj. Volume : 1.000
MW : 774.88

Cat. No.: RP10408
Lot No.: PJ00309510

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp: 250°C
Block Temp: 200°C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN