


MS Spectrum

Positive



Sample Information Item

Date and Time : 2005-4-7 12:57:33  
User : xie  
Sample : RP10447(PEL-13-2)  
Inj. Volume : 1.000  
MW : 1671.98

Cat. No.: RP10447  
Lot No.: PL00215506

Probe: ESI  
Nebulizer Gas Flow: 1.5L/min  
CDL: -20.0v  
CDL Temp: 250°C  
Block Temp: 200°C

Probe bias: +4.5kv  
Detector: 1.5kv  
T. Flow: 0.2ml/min  
B. conc: 50%H2O/50%ACN

m/z